# Deciphering Dysphagia with Ampcare's Effective Swallowing Protocol (ESP™)

Training for the Use of Ampcare's ESP in the Treatment of Dysphagia

**PRESENTED BY:** 



Rick McAdoo, M.S., CCC-SLP Russ Campbell, PT Ronda Polansky, M.S., CCC-SLP





### Presenters

Russell Campbell, PT Ronda Polansky, M.S., CCC-SLP Rick McAdoo, M.S., CCC-SLP

- Russ Campbell practices physical therapy and specializes in the areas of neurology, orthopedics and gerontology; he brings over 30 years of extensive experience in using therapeutic modalities to maximize patient outcomes and is licensed in Texas and Illinois.
- Ronda Polansky is an active clinical Speech Pathologist and co-owner of DiagnosTEX, a mobile modified barium swallow study company in Dallas/Fort Worth, Texas. For over 25 years she has focused her interests on patients with dysphagia.
- Rick McAdoo is recognized for his expertise in working with the adult neurogenic population with specialized emphasis in dysphagia. He was the first to show laryngeal elevation using transcutaneous electrical stimulation on both a computerized laryngeal analyzer and under video fluoroscopy in 1997.
- Collectively Russ, Ronda and Rick are the founders and co-owners of Ampcare, LLC, which was established to develop therapeutic neuromuscular stimulation technologies for the treatment of dysphagia. Among other accomplishments, Ampcare has built a strong patent portfolio, generated clinical data and developed an FDA-cleared therapy system known as the Effective Swallowing Protocol (ESP<sup>™</sup>) for the treatment of dysphagia. Ampcare is committed to delivering positive patient outcomes and innovative products specifically designed to meet the needs of the patient, medical professional and payer by improving the options and availability of dysphagia treatment techniques.

Conflict of Interest Disclosure: As required, included are disclosures by presenters of this continuing education program regarding any significant conflict of interest including financial interest(s) or other affiliations(s) with commercial organizations(s) and any endorsement of products that may have a direct or indirect interest in the subject matter of the program.

This course will focus on Ampcare products and will not include information on other similar or related products. Course offerings will cover information that pertains to the use of the Ampcare-patented methodology. All of the concepts that are presented today are from Ampcare's patent, clinical data, research and development. Ampcare, LLC benefits from the sale of this course and product. Russ Campbell and Rick McAdoo in return have benefited financially from a paid salary through Ampcare. Being co-owners of Ampcare, LLC all three presenters have a financial interest in the products or services described, reviewed, evaluated or compared in this presentation.

Course Financial and In-Kind Support Disclosure: No relevant financial or in-kind support to disclose

The assignment of continuing education hours by approving bodies (ASHA) does not imply endorsement of specific course content, specific products, or clinical procedures by those approving bodies.

Copyright Permission: This material may not be reproduced in any form without the written permission from the program speaker(s)/author(s) listed above.

















6% of Overall Global Population (8B  $\Rightarrow$  480M)

36.5% in Hospital Setting

42.5 % in Rehabilitation Setting

50.2% in Nursing Home Setting

Rivelsrud et al. 2023

mpcare

5













































Patients with dysphagia because of decreased hyolaryngeal excursion and or reduced laryngeal vestibule closure times who receive ESP improved at a higher rate than those patients receiving traditional therapy.

After more than 15,000 therapy visits this protocol is a safe adjunct to indirect treatments for pharyngeal dysphagia with impaired hyolaryngeal excursion and or laryngeal vestibule closure times.

Cochrane Database System Review Conclusions on Swallowing Therapy for Dysphagia

• Authors' conclusions:

npcare

19

• Moderate- and low-quality evidence suggests that swallowing therapy did not have a significant effect on the outcomes of death or dependency/disability, case fatality at the end of the trial, or penetration/aspiration score.

#### • How does this pertain to Electrical Stimulation Research?

- Nearly all of these studies used parameters that are primarily sensory stimulation or TENS. Suiter et al., 2006, Humbert et al., 2006 Bulow et al., 2008
- Some are looking at the immediate effect electrical stimulation might have, which is not rehabilitation. *Clark et al., 2009, Humbert et al., 2015, Ogura et al., 2022*
- Some are using electrode placement that would not follow the rules of NMES and do not address any specific mechanical cause.

mpcare

Ludlow et al., 2007, Bath et al., 2016, Bath et al., 2020





 Ampcare's ESP over only the submental region at rest produced 17% of the extent of hyoid elevation and 4.5% of laryngeal excursion that normally occurs during swallowing.

 Research Utilizing ESP

 Villizing ESP

 Patient at Rest

 Patient Research Produced 17% of the extent of hyoid elevation and 4.5% of laryngeal excursion that normally occurs during swallowing.

 Ludlow et al., 2010





## UK Research Utilizing ESP 2014 Phase II – Pilot Study

- Randomized Controlled Pilot Trial:
  - 30 CVA patients with dysphagia > 1 month
    - Control group (15) usual dysphagia care
    - Intervention group (15) ESP for 30 minutes, 5 times/week for 4 weeks
    - Outcome Measures
      - Functional Oral Intake Scale (FOIS) 57% control, 75% ESP
      - Penetration/Aspiration Scale (PAS) 17% control, 58% ESP for Diet
      - SWAL-QOL (Quality of Life) 42% control, 100% ESP

Sproson et al., 2018

npcare























- Some SLPs have had limited education on the anatomy and physiology specific to a swallow but assess and treat these muscles and cranial nerves.
- Three phases of the swallow work together to produce a functional swallow:
  - Oral Preparatory/Transport
  - Pharyngeal
  - Esophageal
- Any breakdown in a phase, muscle, or cranial nerve can cause dysphagia.





























A-P propulsion/transit of the bolus

Normal movement of

the tongue is essential

for carrying out the

tasks of the oral phase

of a swallow.

The tongue provides

anchoring during hyoid

bone elevation and UES

opening. It connects to

the hyoid bone,

pharynx, and epiglottis.

Begins when the tongue contracts against the hard palate.

mpcare

39













Swallow Problem	Muscle Group	<b>Cranial Nerve</b>
Acceptance	Orbicularis Oris, Zygomaticus	CN VII
Anterior Loss	Orbicularis Oris	CN VII
Mastication	Masseter, Pterygoids, and Temporalis	CN V
Oral Pocketing or Scattering	Buccinator and intrinsic and extrinsic lingual muscles	CN VII, XII
Anterior-posterior propulsion	Intrinsic and extrinsic lingual muscles – Genioglossus, Styloglossus	CN XII
Poor lingual pressure to drive the bolus through the PES/UES	Intrinsic and extrinsic lingual muscles – Hyoglossus	CN XII

























## Infrahyoids (strap muscles) HYOLARYNGEAL DEPRESSORS

Omohyoid – depresses hyoid

Sternohyoid – depresses hyoid

Sternothyroid – depresses thyroid

Thyrohyoid – shortens the distance between the thyroid and hyoid bone











Swallowing Problem	<u>Muscle Group</u>	Cranial Nerve	
Poor velopharyngeal seal	Tensor Veli Palatini, Pharyngeal Constrictors, Levator Veli Palatini, Musculus Uvulae	CN V, IX, X, XI	
Decreased closure of the larynx	Aryepiglottic, Lat. Cricoarytenoid and Transverse Arytenoid	CN X, XI	
Weak pharyngeal constrictor contraction	Superior, Middle, and Inferior Pharyngeal Constrictors	CN IX, X, XI	
Decreased anterior superior elevation of the hyolaryngeal complex	Anterior Digastric and Mylohyoid	CN V	
Failure of opening the UES/PES	Inferior Pharyngeal Constrictor, Cricopharyngeus, Superior Longitudinal Esophageal muscle	CN IX, X, XI	










































IVIUSCIE FIDER TYPES Muscle fibers can be grouped according to the kind of tissue they are found in:		
Muscle Type	Location	Contraction Contro
Cardiac	In the walls of the heart	Involuntary
Smooth	In the walls of hollow internal structures (e.g., esophagus)	Usually Involuntary
Skeletal	At the attachment to bones	Voluntary















Fiber Type	Туре І	Type IIa	Type llb
Contraction time	Slow	Fast	Very Fast
Endurance	High	Intermediate	Low
Fatigue	Slow	Intermediate	Fast
Power	Low	High	Very High
Motor neuron size	Small, Narrow	Large	Very Large
Function	Static, Postural	Dynamic	Dynamic, Explosive
Energy Source	Oxygen	Oxygen, Glycogen	Glycogen
Color	Red	Pink	White















































### Electrical Stimulation for Muscle Re-education

Process by which electrical stimulation works has been well documented and studied.

Extrinsic muscles of the swallowing mechanism bear no histological difference to any other skeletal muscles typically treated by PT and OT.

mpcare

91

	Prevent	Prevent disuse atrophy	
	Increase	Increase range of motion	
ladiosticas	Re-educate	Re-educate muscle function	
indications	Decrease	Temporarily decrease spasticity	
	Serve	Serve as an electrical orthosis	
	Increase	Increase local blood circulation	
mpcare	L	l]	





























































### Ampcare's ESP Effective Swallowing Protocol

- Symmetrical Biphasic waveform.
- Frequency of 30 pulses per second.
- **Phase Duration** of 50 microseconds but if the patient can tolerate maximum intensity we recommend 250 microseconds.
- Ramp up of 1 second; Ramp down of 0 seconds.
- **Duty cycle** of 5 seconds on, 25 seconds off initially, but progress to 5 seconds on, 20 seconds off and finally 5 seconds on and 15 seconds off.
- Treatment time for 30 minutes.
- Intensity or amplitude to tolerance, but preferably high enough for a sustained motor contraction.
- Patient to perform indirect exercise while electrical stimulation is on.



111





## Geniohyoid muscle contraction during stim



Device: ITO postim Mode: NMES Frequency: 30pps Phase duration: 50µs Ramp up time: 1s Stim control: Foot switch





Dr. Tatsuyuki Fukuoka (Hiroshima International University)

113

























Ampcare

119





# E Series Electrode Use

- 1. Wash the treatment area with soap and water, rinse, and dry well.
- 2. Skin prep wipe may also be used. Allow the skin to dry for at least 45 seconds.
- 3. Remove electrodes from the liner by lifting a corner and pulling gently.
- 4. Place electrodes on the correct body part to be stimulated.
- 5. A strap or **Restorative Posture Device (RPD)** can be used to hold electrodes in place to increase conduction.
- 6. Attach electrodes to the lead wires.
- 7. Make sure the unit is off and attach the lead wire to the stimulator.
- 8. Turn on the stimulator and increase the intensity to tolerance but preferably high enough for a tetanic contraction.























### Modifications to the RPD

• Proper fit will depend on preexisting conditions, comfort and therapeutic goals.



127













Program	ESP1	ESP2	Manual 3
Frequency (Hz)	<b>30 (</b> †	fixed)	5-50 adjustable in 1 Hz increments
Phase Duration (µsec)	50 (fixed)	250 (fixed)	50-250 adjustable in 50 μsec increments
Ramp Up (sec)	1 (fixed)		1 (fixed)
Ramp Down (sec)	0 (fixed)		0 (fixed)
Cycle ON Time (sec)	5 (fixed)		
Cycle OFF Time (sec)	Selectable values of 25, 20, or 15		
Program Duration (min)	30 (fixed)		5 – 30 adjustable in 5 min increments
ntensity	0 – 20 adjustable in increments of 1*		adjustable in increments of 1*





133







Example of oral motor impairments:

- Facial drooping
- Drooling
- Incomplete mastication
- Oral pocketing
- Swallowing and speaking

This can be due to facial weakness or abnormal muscle tone.











Facial Nerve (VII)	
Muscle	Function
Orbicularis Oris	Controls lip closure, pursing movements – intake through a straw or kiss someone. Allows you to say "puh", "MMM", /w/ or "ooo" sounds.
Buccinator	Compresses the cheeks against the teeth to prevent oral pocketing of food, also used in the act of blowing air between the lips as in whistling or to increase oral pressure, speech intelligibility and volume.
Zygomaticus Major	Draws the angle of the mouth upward and outward as in smiling and for oral acceptance. Allows you to say "eee" sounds.







Muscle	Function
Masseter	Provides jaw closure and range of motion to masticate solids and assist with articulation.
Anterior Digastric, Mylohyoid	Hyolaryngeal elevation and protraction and mara affect voice pitch. Assists with mandibular depression


























# **Treatment Guidelines**

Voluntary exercise should occur simultaneously with the electrical stimulation (e.g., Resistive Chin-to-Chest, Mendelsohn maneuver, Effortful swallow, and oral motor exercises)

Maximal evoked contractions as tolerated to facilitate explosive strength (type II muscle fibers)

As patient fatigues attempt to modify the activity to continue to evoke contractions



147











# **Daily Treatment Note**

### Example:

Patient received ESP1 to the suprahyoid musculature x 30 minutes, 5 sec on/(25) sec off duty cycle with intensity <u>setting</u> of <u>11</u> while performing 10 minutes each of chin-to-chest with effortful swallows (18/20), Mendelsohn (16/20), and jaw open/jaw closed with effortful swallows (14/20).



151

# <section-header><text><text><text><text><text><text><text>















# Post Test

- 2. The indications for NMES include the following EXCEPT:
  - a) Prevent disuse atrophy
  - b) Increase range of motion
  - c) Re-educate muscle function
  - d) Decrease blood circulation



157









LAB 011617 Rev. F











- 8. Electrical stimulation applied to the skin at these levels will penetrate deeper and depolarize nerve endings.
  - a) Increased intensity, increased phase duration.
  - b) Decreased intensity, decreased phase duration.
  - c) Intensity or phase duration do not factor in the depth of penetration of the current.
  - d) Works best when applied to cartilage.































# Deciphering Dysphagia



# Deciphering Dysphagia





Permission for AMPCARE, LLC to reproduce these illustrations from <u>Grays Anatomy, Vol. 36</u>, 1980, Williams et al granted by Elsevier Ltd, The Netherlands.





Permission to reproduce the illustrations from Dr. Wesley Norman was granted to AMPCARE, LLC, Copyright© 1999 by Wesley Norman, PhD, DSc



# **Strength Duration Curve**

### Pacemaker/Defibrillator Protocol

It is always best practice to notify the patient's cardiologist or cardiac clinic (if the cardiologist is no longer seeing this patient) on a patient by patient basis, to clarify if using NMES would be appropriate. Depending on the type of cardiac device the patient may or may not be appropriate for NMES.

If it is approved by the physician you should consider taking the appropriate steps as discussed in the CE Training. Certain pacemaker companies manufacture specific types of pacemakers/defibrillators that can utilize an external magnet to lock the pacemaker signal in place and temporarily prevent the defibrillator from firing. You will need to ask the physician to provide you with the serial number of the device, the company representative's name and phone number in your area to set up an appointment, if the physician agrees on the NMES plan of care on a patient with a pacemaker. The pacemaker company representative may monitor the patient with portable ECG/telemetry during treatment, and instruct you on the process of using the magnet with ESP. If this process is completed and there is no change in the ECG, the company representative may give you a magnet for continued use. It may be that the the company representative just provides you with a magnet without this ECG/Telemetry step. For a step by step process of the Pacemaker/Defibrillator protocol, see below.

### Pacemaker/Defibrillator Protocol:

- 1. Contact Cardiologist or Physician at the cardiac clinic that cares for the patient:
  - Explain that you are treating his/her patient for dysphagia and would like to use NMES.
  - Discuss placement of electrodes over the chin in the submental region using a frequency of 30 Hz and phase duration of 50 μsec to 250 μsec.
  - Request that you would like to set up an appointment with the pacemaker representative to see if an external magnet would be appropriate to prevent electromagnetic inference (EMI).
  - If physician indicates that this is unnecessary, ask for documentation, such as a prescription that indicates "Appropriate to proceed with NMES."
- 2. Contact Pacemaker Company
  - Call number on back of patient's pacemaker card and ask to speak with a representative in your area.
  - Request an appointment for electromagnetic interference (EMI) testing.
  - Give requested information, such as reason, location, etc.
- 3. Discussion with Pacemaker Representative
  - Explain that you want to use NMES with the patient and that you have spoken with cardiologist.
  - Schedule time with representative to allow both of you to be present with the patient.
  - The EMI testing should be a free service.
- 4. Schedule the Appointment
  - Company representative will set the patient up to test the pacemaker/defibrillator for EMI.
  - Place the electrodes submentally and be ready to begin the treatment.
  - The representative will place a magnet over the pacemaker/defibrillator to prevent the defibrillator from firing, but still allows the pacemaker to function.
  - Explain to the patient that you need to increase the intensity as high as you can to ensure there will not be a problem during later treatments.
  - Set the device with standard parameters indicated for treating that particular patient.
  - Turn on and progress the intensity as high as they can tolerate; the representative will have a digital display to tell if there is any EMI; if so discontinue and do not proceed with use. If not, you can proceed.
- 5. Documentation
  - Ask representative for a printout of the ECG.
  - Sign it, along with signature from the representative.
  - Write the date and the parameters used and document in the chart as well
  - Example: "Medtronic representative John Smith provided pacemaker/defibrillator EMI testing for NMES in the submental region for potential treatment of dysphagia. Equipment tested at frequency of 30 Hz, phase duration of 50 µsec, intensity of 55 milliamps with placement in submental region. No EMI noted."
  - Send a copy of findings to cardiologist or physician at the cardiac clinic that cares for the patient.

# Treatment Algorithm for ESP™



# **Types of Electrical Stimulation**

	Direct Current	lonto- phoresis	Edema, Pain	10Hz- 50kHz	Continuous or pulsed	1-2	Bipolar	Mono- or biphasic, or quadratic	1-4mA
nts		TENS	Pain	1-120Hz adj	50-450µs	1-2	Bipolar or Quadripolar	Modified Asymmetrical Biphasic Square Wave	1mA-70mA into 500Ω
ectrotherapy Currer	<b>Pulsed Current</b>	Neuromuscular Electrical Stimulation (NMES)	Prevent disuse atrophy, increase ROM muscle re-education	1-160Hz adj	<b>50-250μs</b>	1-4	Bipolar or Quadripolar	Asymmetrical or Symmetrical Biphasic Square Wave	1-100mA into 500Ω
Ш		НVG	Edema, Wound	Care 1-120Hz adj	100µs	1-2	Bipolar	Mono- phasic Dual Spiked	Wave 1-300V
	Alternating Current	IFC Premod Russian	Pain Edema	Carrier Freq 4000Hz, Modulated Freq 4001- 4160Hz	125µs	1-4	Bi or Quadripolar	Symmetrical balanced sine	1-60mA into 500 <u>0</u>
			Indication	Frequency	Phase Duration	Channels	Set-up	Waveform	Intensity

# Deciphering Dysphagia with E-Stim

Neuromuscular electrical stimulation can help you successfully treat dysphagia by **jump-starting** the biomechanical movement of swallowing.

> By Russ Campbell, LPT, Ronda Polansky, MS, CCC-SLP, and Rick McAdoo, MS, CCC-SLP

EATING AND DRINKING COMES NATURALLY FOR MOST OF US. For people with dysphagia, swallowing a bite or a sip is difficult—and in some cases deadly.

This complex condition is not easily solved. But a new treatment approach using neuromuscular electrical stimulation (NMES) may hold a hint to the riddle of dysphagia. NMES has proven successful in treating a variety of neu-

romuscular and musculoskeletal problems, provided the peripheral nervous system is intact. It has been reported to increase strength and range of motion, facilitate weak contractions due to upper-motor neuron lesions or disuse atrophy, and re-educate muscles.<sup>1-5</sup>

In dysphagia rehabilitation, NMES can jumpstart the biomechanical movement of the anatomical structures associated with swallowing. These include the mylohyoid, geniohyoid and anterior belly of the digastric muscles, which are primarily responsible for anterior and superior movement of the hyoid bone during a

Digastric muscle

Anterior Belly ..

Hyoid bone

swallow. This movement of the hyoid and laryngeal elevation is vital in preventing aspiration. The superior movement of the larynx helps bring the airway safely away from the path of the bolus.

A preliminary study performed at our facility has located a motor point that can be stimulated to promote laryngeal elevation, which is often a problem in dysphagia. A physical



therapist set the parameters for the electrical stimulation, and a speech-language pathologist analyzed the laryngeal elevation.

A Computerized Laryngeal Analyzer (CLA),<sup>6</sup> approved by the FDA for diagnosing and treating swallowing disorders, was used as the measuring instrument.

The technology is noninvasive; it measures the biomechanical movements during the swallow and automatically calculates and displays the duration of the laryngeal motion. This allows the clinician to record the swallow pattern, intensity and timing and assign them a quantitative value. *Continued on page 67* 

Mylohyoid muscle

In dysphagia rehabilitation, neuromuscular electrical stimulation can jump-start the biomechanical movement of the anatomical structures associated with swallowing.

### DYSPHAGIA

### DYSPHAGIA

Continued from page 64

The CLA uses a thin film piezoelectric transducer, software and a high-speed computer to provide on-line and real-time display. The graph displays superior and anterior planes of movement upward and inferior, and posterior planes of movement downward. Wave forms on the CLA correspond to the laryngeal muscles' activity during the swallow.

To use NMES to promote laryngeal elevation while recording it using the CLA, clinicians attach a piezoelectric transducer to the subject's thyroid cartilage. Three dry swallows are monitored by the CLA.

Next, two moveable electrodes connected to a neuromuscular electrical stimulator are placed over the region of each belly of the anterior digastric and mylohyoid muscles. The benefits of targeting these muscles are twofold: They voluntarily assist the end of the oral stage through the elevation of the tongue, and the muscles contract during the pharyngeal stage, assisting elevation of the larynx.<sup>7</sup>

The electrodes are placed on an area between the oral cavity and the circumference of the neck to decrease the risk of adverse effects as listed by the FDA. The electrodes are housed in plastic caps with a moist sponge of water conducting the current to the muscles.

The patient is asked to hold the NMES unit while pressing the constant stimulation (CS) button and increasing the intensity until a vibration is felt. The physical therapist moves the electrodes on the muscle until he finds a motor point and observes visual tetanizing contraction. A motor point is where a motor nerve enters the muscle it innervates.<sup>8</sup>

To achieve a visible contraction, the patient may need to increase the intensity. If the intensity becomes uncomfortable without visible contraction, the electrodes may be over a sensory nerve point. Once the motor points are found, the patient will increase the intensity to maximum tolerance and release the CS button.

The physical therapist then takes the NMES unit from the subject and stimulates the area three times for approximately four seconds each, while the patient holds the electrodes and maintains a neutral head and neck position. Each subject is instructed to look straight ahead and not assist or resist the contraction in any way during each



four-second trial.

Approximately 15 to 30 seconds of rest time are given between stimulation trials while the speech-language pathologist assesses the elevation in millivolts on the CLA. The parameters include:

• a pulse rate of 30 Hz with a preset pulse width to produce a tetanizing muscle contraction with minimal muscle fatigue<sup>5</sup>

• a .4-second ramp to allow for rise in intensity

symmetrical biphasic waveform to efficiently stimulate both electrodes.<sup>9</sup>

After the .4-second ramp of stimulation, an initial recruitment of the platysma creates a downward motion of the larynx. The platysma—a broad thin layer of muscle situated on each side of the neck immediately under the superficial fascia—draws the lower lip and corner of the mouth to the side and down. When moved forcefully, it expands the neck and draws the skin upward.<sup>10</sup> After involving the platysma, the suprahyoidal musculature then elevates the larynx. This pressure is maintained with some fasciculations until the stimulation is released at the four-second frame on the CLA graph.

Properly placing the two electrodes over the motor points of each belly of the anterior digastric and/or mylohyoid muscles elicited mea-▷

# The Software Solution That Integrates

Bar Code Scanning Electronic Billing Utilization Breakdowns Multi Doctor/Multi Clinic



We offer our customers total business solutions. Our wireless hand held scanner enables you to enter information as you go from room to room or clinic to clinic. Info can then be downloaded to update the ledger and prepare for automatic billing.

Raintree Systems Inc.

"The Software Solution You'll Never Outgrow" Raintree now integrates with TherAssist

### For Dealer Nearest You Call 1-800-333-1033

Circle 43 on Info Card November 1998 • advance for Directors in Rehabilitation 67

# **Busy day?** Few minutes to spare?

## Make them count.

Subscribe to ADVANCE online at advanceweb.com Or call 1-800-355-1088.



# How do you view your outpatients?



### The picture becomes CLEAR with LIFEware. **Using the Profile Report lets you:**

• Reduce Clinician Assessment Time • Focus on Areas of Deficit • Standardize Among all Clinicians • Improve Treatment Quality Document Progress to Payers

Using LIFEware<sup>®</sup> to track outpatient outcomes can do more than simply aggregate data so you can see how you're doing vs. other facilities. By giving a comprehensive picture of how the patients see themselves, LIFEware<sup>®</sup> can help you formulate a more effective treatment plan that in turn, should be more cost-effective.

### Check Out the LIFEware Online Demo by Visiting Our Website: www.lifeware.org

LIFEware<sup>®</sup> • UDSmr • The FIM System Phone: 716-829-2076 ext 32 • E-mail: info@lifeware.org

### DYSPHAGIA

surable movement of the larynx.

The movement recorded averaged a mean of 40.04 mV on 25 people through using NMES only. The participating subjects recorded a mean of 52.6 mV on a normal dry swallow, with consistent waveforms. These two results denote that NMES elicited 76 percent of the mean laryngeal elevation of the normal swallows.

The initiation of movement achieved through electrical stimulation recorded on the CLA resembles a pattern that demonstrates laryngeal elevation while maintaining stimulation for four seconds.

Movement of the larynx was displayed when the stimulation began. During the stimulation, the larynx is held in position; therefore, there are no changes in pressure from the thyroid cartilage and the piezoelectric transducer. This is demonstrated by the return of the graphic display to the 0 millivolt level.

At the stimulation release, a new change in pressure occurs when the larynx returns from the elevated position. It is significant that a sustained elevation was achieved by stimulation only (no swallow initiation) and required no understanding or mastering of a special instruction or technique. Patterns and millivolts were comparable to each person's normal laryngeal elevation pattern and the amount of intensity tolerated (intensity range 15-25).

These preliminary findings indicate a positive breakthrough in the

### Electrical stimulation can promote laryngeal elevation without a volitional

response from the subjects . . . One day, this research could *yield the final* clue to deciphering dysphagia.



use of NMES in dysphagia rehabilitation. Electrical stimulation can promote laryngeal elevation without a volitional response from the subjects. In the future, this process could be used with other instrumental examinations to investigate epiglottic closure and esphageal sphincter elongation. One day, this research could yield the final clue to deciphering dysphagia.

### References

1. Takebe, K., & Basmajian, J.V. (1976). Gait analysis in stroke patients to assess treatment of foot drop. Archives of Physical Medicine and Rehabilitation, 57:1, 305-310.

2. Lee, K.H., & Johnston, R. (1976). Electrically induced flexion reflex in gait training of hemiplegic patients: Induction of the reflex. Archives of Physical Medicine and Rehabilitation Medicine, 57:1, 311-314.

3. Cranston, B. et al. (1977). Improvement of gait following functional electrical stimulation: Investigations on changes in voluntary strength and proprioceptive reflexes. Scandinavian Journal of Rehabilitation Medicine, 9(1), 7-13.

4. Merletti, R. et al. (1979). Clinical experience of electronic peroneal stimulators in 50 hemiplegic patients. Scandinavian Journal of Rehabilitation Medicine, 11(3), 111-121

5. Benton, L., Baker, L. et al. (1981). Functional electrical stimulation. A practical guide, 2nd Edition. Ranchos Los Amigos, Rehab Engineering Center.

6. H.I.L. Research, Ltd. (1995). CLA: Computerized laryngeal analyzer. Tel Aviv, Israel.

7. Kendall, F., & McCreary, E. (1983). Muscles testing and function, 3rd Edition. 266. Baltimore, MD: Williams & Wilkins.

8. Mosby's Medical and Nursing Dictionary. (1983). St. Louis, MO: Mosby.

9. Respond Select (1996). NMES. St. Paul, MN: EMPI, Inc.

10. Webster's Medical Desk Dictionary. (1986). Springfield, MA: Merriam Webster Inc.

Russ Campbell, LPT, is senior physical therapist at Bowie Rehabilitation Services and Ronda Polansky, MS, CCC-SLP, is senior speech pathologist; Rick McAdoo is operations manager at Bowie Speech and Hearing Inc., Fort Worth, Texas.

68 NOVEMBER 1998 • advance FOR DIRECTORS IN REHABILITATION LAB 011617 Rev. F

Circle 44 on Info Card

Neuromuscular e-stim offers another option to treat decreased laryngeal elevation.

### **BY RONDA POLANSKY, MS,** CCC-SLP, RUSS CAMPBELL, PT, & RICK MCADOO, MS, CCC-SLP

ith age comes a range of problems, from diminished sight and poor hearing to altered mobility and balance difficulties. In addition to those issues, as many as 70 percent to 90 percent of elderly patients have some degree of swallowing dysfunction.1

Dysphagia affects 15 million people in the United States and 45 percent to 87 percent of residents in long-term care.2-4 Aspiration pneumonia is the fifth leading cause of death in the United States.<sup>5</sup> The cost to treat dysphagia that results in pneumonia is estimated to exceed \$3 billion each year. Approximately 150,000 nursing home patients require hospitalization for pneumonia each year.1

The amount of time a speech-language pathologist spends caring for patients with dysphagia has increased over the past 20 years due to more referrals and better awareness among health professionals. Speech-language pathologists should play a critical role during the dysphagia recovery process.

### **DEALING WITH DYSPHAGIA**

The symptoms of dysphagia vary with its cause and can have a devastating effect on quality of life. Direct and indirect treatments include oral-motor exercises, hvolarvngeal and Shaker exercises, the Mendelsohn maneuver and thermal stimulation.

These treatments are often difficult to perform with the geriatric population when there is a comorbidity with Alzheimer's disease, stroke

or Parkinson's disease. Decreased laryngeal elevation is a common functional abnormality in pharyngeal dysphagia. Laryngeal elevation is important in airway protection to prevent penetration and aspiration.

Electrical stimulation (e-stim) has been used for years as a physical therapy modality. Although transcutaneous e-stim has received more attention as an adjunct to dysphagia treatment recently, clinicians don't know much about its effect on swallowing physiology. E-stim assists hyolaryngeal elevation and may increase sensory input to the central nervous system to enhance swallowing.6

Using neuromuscular electrical stimulation (NMES) to treat laryngeal elevation provides speech-language pathologists with another treatment option to help patients improve swallowing disorders, and enable them to return to feeding and a regular diet. When NMES is applied to the skin at low current levels, it activates sensory nerve endings in the surface layers, which provides sensory feedback to the peripheral and central nervous systems.

With increased intensity or pulse width, the electric field penetrates deeper and depolarizes nerve endings in muscles lying beneath the skin surface. The field also spreads with diminishing density to produce a muscle contraction, provided the peripheral nervous system is intact.

The U.S. Food and Drug Administration has cleared NMES for use for muscle re-education, to prevent or retard disuse atrophy and relax muscle spasms. In addition, NMES is FDAapproved to increase local blood circulation, maintain or increase range of motion and for immediate post-surgical stimulation of calf muscles to prevent deep

vein thrombosis.

Using e-stim for this program is based on a standard muscle re- pathologists education protocol for small muscle groups. It's also applied to pre- critical role vent or retard disuse during the atrophy of the suprahyoid musculature due dysphagia to inactivity.

To treat diminished hyolaryngeal elevation, process. electrodes are placed

Speechlanguage should play a recovery

www.advanceweb.com/rehab

LAB 011617 Rev. F

ADVANCE FOR DIRECTORS IN REHABILITATION . MAY 2008 39

# **Ora-Light**



Ora-Light provides a balanced program of lip and tongue excercises using an innovative set of practical tools. Ora-Light helps to maximize rehabilitation outcomes of patients with oral motor deficits including dysarthia, dysphagia, apraxia, articulation disorders, oral cancer resections and laryngectomee rehabilitation. For a free demonstration DVD call us today at 207-741-2443.

**Speech Pathology Associates, LLC** P.O. Box 2289 South Portland, ME 04116

www.chewytubes.com

### DYSPHAGIA

on the submental area of the submandibular triangle, rather than the contraindicated portion of the anterior neck area, such as the carotid sinus or laryngeal region. This location is safer and more physiologically plausible.

By using standard equipment for this purpose, clinicians can modify the parameters to meet a patient's individual needs. This can be done with adequate electrotherapy training.

Used properly, NMES is safe and effective to treat poor hyolaryngeal elevation. Through continued research, this approach has been proven under videofluoroscopy, as well as fiberoptic endoscopic evaluation of swallowing (FEES) supporting laryngeal elevation, tongue base retraction and improved by-mouth intake.<sup>7</sup>

### **PUTTING IT TO THE TEST**

A clinical trial evaluated the efficacy of using e-stim to the anterior digastrics and mylohyoid musculature (submandibular region) to promote improved laryngeal elevation and improved swallow function and safety. The rationale was that range of motion of the larynx during swallowing affects airway protection, and anterior motion of the hyolaryngeal complex is essential to successful function of the cricopharyngeal muscle. The trial focused on using these e-stim parameters with a restorative postural device to promote optimal posture and electrode conduction, combined with electrodes to re-educate and strengthen muscles related to laryngeal elevation.

The trial involved patients at long-term care facilities in Texas. Patients in group 1 received at least 20 therapy days of the NMES protocol, along with traditional therapy (79 patients). Patients in group 2 received traditional dysphagia therapy, which would have been per-



formed prior to NMES modality training or for patients who refused the NMES therapy protocol (46 patients).

Diagnoses included cerebrovascular accidents, Parkinson's, Alzheimer's, dementia, pharyngeal cancers and pneumonia from a pharyngeal dysphagia.

To date, 125 patients evaluated by modified barium swallow study (MBSS) had impaired laryngeal elevation as a primary or secondary dysfunction, causing aspiration or risk of aspiration to the degree that diet changes were necessary. A swallow severity scale was established to determine the diet.

The subset of patients who could tolerate at least 20 days of traditional dysphagia therapy while using the NMES protocol were included in the analysis as successfully completing the protocol. They were compared with 46 patients who received only traditional dysphagia therapy, but whose chart reviews noted they exhibited dysphagia with decreased laryngeal elevation as diagnosed from MBSS under fluoroscopy.

The charts were evaluated as to the number of patients who had an improved swallowing severity scale. The six scale levels are: nothing by mouth (0); therapeutic intake only (1); pleasure feedings only (2); modified diet including thickened liquids and puree or mechanical soft with strategies (3); strategies only, no alternate method of intake (4); and normal swallowing function (5).

The scale level improved from 2.33 to 3.72 in the NMES subgroup that received at least 20 days of the NMES protocol, and from 2.52 to 2.6 in the traditional therapy group. Not all of the patients were able to achieve a period of at least 20 days in the traditional therapy group.

The average number of therapy days was 36.79 in the NMES subgroup and 19 in the traditional study group. Three-fourths of the patients who received at least 20 days of the NMES protocol had a diet upgrade.

The results of this clinical trial suggest that patients who present with dysphagia, due in part to diminished laryngeal elevation, and receive NMES to the laryngeal elevators improve diet upgrades and swallow function at a higher percentage. Throughout the trial, there were more than 5,814 therapy visits using this protocol. This suggests that it's a safe adjunct to treat patients with pharyngeal dysphagia.

As research expands in this area, clinicians will be able to determine more about the physiological effects on the neck muscles and swallowing functions.

### References

1. Murray, J., & Carrau, R. (2006). Clinical management of swal-

www.advanceweb.com/rehab

lowing disorders. San Diego, CA: Plural Publishing. 2. Logemann, J. (1991). Manual for the videofluoroscopy study of

swallowing. Austin, TX: Pro-Ed. 3. Langmore, S. (1994). Predictors of aspiration pneumonia.

Dysphagia, 9, 90-95. 4. Steele, C.M., Greenwood, C., Ens, I., et al. (1997). Mealtime

difficulties in a home for the aged: Not just dysphagia. Dysphagia, 12, 43-50.5. Langmore, S.E. (1991). Managing the complications of

aspiration in dysphagic adults. New York: Thieme Medical Publishers.

6. Campbell, R., Polansky, R., & McAdoo, R. (1998). Deci-

**© GRAVITY** 

phering dysphagia with e-stim. ADVANCE for Directors in Rehabilitation, 7(11), 64, 67, 68.

7. Campbell, R., Polansky, R., McAdoo, R. (1998). Promoting laryngeal elevation with e-stim. ADVANCE for Speech-Language Pathologists & Audiologists, 8(25).

Ronda Polansky, MS, CCC-SLP, Russ Campbell, PT, and Rick McAdoo, MS, CCC-SLP, are partners and co-owners of Ampcare LLC in Fort Worth, Texas. They can be reached at (682) 561-2444 or ampcare@juno.com



Full implementation training. Included.

The Cash Practice Paradigm

ringer hand and manufacture hand the

See us at APTA PT San Antonio, June 11-14 Booth 724 See us at NATA St Louis, June 18-20 Booth 1337

Serving the rehabilitation community with Total Gym since 1985.

www.ampcarellc.com



Card

on FREE Information

33

Circle

800 541 4900 Call to receive your Free Presentation DVD. Visit our website at efisportsmedicine.com

ADVANCE FOR DIRECTORS IN REHABILITATION • MAY 2008 41

An investigation into the ef	fects of neuro-muscular electrical s in patients presenting with persiste	timulation (NMES) on swallowing ent dysphagia.
Speech and Langu Sc	Dr Sue Pownall and Professor Pam Ender age Therapy Department Sheffield Teaching Hospital NH hool of Health and Related Research, University of Sheffi	rby S Foundation Trust. Correspondence to: e <i>ld.</i> sue.pownall@sth.nhs.uk
1) Background	3) Results	4) Results
Decreased laryngeal elevation is common among patients with swallowing difficulties. It has been hypothesised that electrical stimulation can assist hyo-laryngeal elevation (Freed et al 2001; Leelamanit et al 2002). Stimulation of laryngeal elevation musculature via NMES in combination with traditional treatments is	5 patients were recruited. Two with diagnosis of stroke, two head and neck cancer and one skull base osteomyelitis. Time of dysphagia ranged from 3 months to two years. All participants completed 20 sessions, described the intervention as acceptable and would strondly recommend it to others.	Patient description of daily oral intake (All had PEG tubes for nutrition) Patient 1 Pre-5 teaspoons fluid or 5 teaspoons soft fork mashed diet Follow-up- Eating full oral intake. PEG to be removed
suggested to re-educate muscles when suffering from reduced laryngeal elevation. The rationale for the technique and parameters follow the rules of electrophysiology for stimulating small muscle groups. Using NMES to treat decreased laryngeal elevation provides SLTs with another treatment option.	All 5 participants showed improvements in their swallowing to varying degrees. All 5 participants showed improvements in their swallowing to varying degrees. Two returned to full oral diet having been on non-oral feeding for 5 months and over 24 months respectively, although one continued to show aspiration with some intake. Both ceased nutritional intake via PEG reducing costs to the NHS.	<ul> <li>Patient 2 Pre- Half pint soup taken in small amounts throughout day Follow-up- As above plus 4x half cup coffee per day. Introducing small amount soft fork mashable foods</li> <li>Patient 3 Pre- small amounts normal diet and fluids Follow-up- Full oral intake. Flushing PEG only to keep patent</li> </ul>
	Three patients introduced an increased amount and range of food consistencies to their oral intake, although they continue to receive their main nutrition via PEG.	Patient 4         Pre- 2-4 cups tea per day           Follow-up- intermittently eating minimal amounts soft fork mashed diet. First solid food in 3 years
<ul> <li>2) Method</li> <li>Feasibility study using a case series design which included 5 patients with stable and persistent dysphagia who presented with</li> </ul>	Functional Oral Intake Scale	Patient 5 Pre- Water only sucked from mouth sponges Follow-up- Sips water, plus small amount soft fork mashed diet
decreased laryngeal elevation diagnosed by Videofluoroscopy (VF). Participants received 4 weeks of daily intervention (5 days per		5) Conclusion
week for 30 minutes). Participant's swallowing was re-assessed by VF and clinical assessment post intervention. A follow-up assessment was completed two weeks after completion of the intervention.	2 Follow-up	Results suggest this technique may be an effective intervention for some patients presenting with persistent dysphagia.
Data was collected using Rosenbek Penetration/aspiration scale, Waxman dysphagia severity scale and Functional Oral Intake Scale. Participants also completed a questionnaire about accentability of the intervention		A pilot study is now planned to investigate the technique using a randomised design.
Data was analysed using descriptive statistics.	<ul> <li>Level 1: Nothing by mouth.</li> <li>Level 2: Tube dependent with minimal attempts of food or liquid.</li> <li>Level 3: Tube dependent with consistent oral intake of food or liquid.</li> <li>Level 4: Total oral diet of a single consistency.</li> <li>Level 5: Total oral diet with multiple consistencies, but requiring special preparation or compensations.</li> <li>Level 5: Total oral diet with multiple consistencies without special preparation, but with special diet with multiple consistencies without special preparation, but with special or diet with multiple consistencies without special preparation, but with special or diet with no restrictions.</li> </ul>	References Freed M.L. Freed L. Chatburn R.L. Christian M. Electriscia stimulation for swallowing disorders caused by stroke. Respir Care 46 (5):466-474, 2001 Leelamanti V. Limsakul C. Geater A. Synchronized electrical stimulation in treating pharyngeal dysphagia. Laryngoscope 112(12): 2204-2210, 2002

The



Archives of Physical Medicine and Rehabilitation

journal homepage: www.archives-pmr.org Archives of Physical Medicine and Rehabilitation 2013;94:2542-8

### **ORIGINAL ARTICLE**

# Measurement of Hyolaryngeal Muscle Activation Using Surface Electromyography for Comparison of Two Rehabilitative Dysphagia Exercises



Christopher R. Watts, PhD

From the Department of Communication Sciences & Disorders, Texas Christian University, Fort Worth, TX.

### Abstract

**Objective:** To investigate the effects of a resistance-based chin-to-chest (CtC) exercise on measures of hyolaryngeal muscle activation compared with a head-lift exercise.

Design: Within-subject, repeated-measures design.

Setting: Academic research laboratory.

**Participants:** Healthy young women (N=20) without a history of dysphagia, cervical spine conditions, neurologic disease, or head/neck cancer (mean age, 22.5y).

**Interventions:** All participants performed an isometric jaw-opening exercise against resistance (CtC) and an isometric head-lift exercise, both targeting activation in the hyolaryngeal (suprahyoid) muscles. The CtC exercise required jaw opening into a chin brace secured against the upper torso for a duration of 10 seconds. The isometric head-lift exercise required lifting and holding the head from a supine position for 10 seconds. The degree to which each exercise activated the suprahyoid muscles was measured using surface electromyography (sEMG).

**Main Outcome Measures:** Microvolts as measured from sEMG sensors placed on the skin surface above the hyolaryngeal muscles (surface of skin above geniohyoid, mylohyoid, and anterior digastric). Dependent variables included the peak microvolts during 10 seconds of sustained contraction and the difference in microvolts from rest to peak contraction for each exercise.

**Results:** Activation in the hydraryngeal musculature as measured via sEMG was significantly greater when participants performed the CtC exercise compared with the head-lift exercise. Measures of peak microvolts during contraction were significantly greater for CtC (t=10.72, P<.001) compared with the head-lift exercise, and difference measures in microvolts calculated between rest and contraction for each exercise revealed a 2-fold increase in hydraryngeal muscular activation for CtC (t=8.27, P<.001).

**Conclusions:** The isometric CtC exercise resulted in greater activation of the hyolaryngeal muscles compared with an isometric head-lift exercise. Results support the need for further investigations to determine whether the CtC exercise has a positive effect as a rehabilitative exercise for clinical populations with dysphagia secondary to upper esophageal sphincter dysfunction where hyolaryngeal excursion is a physiological impairment.

Archives of Physical Medicine and Rehabilitation 2013;94:2542-8

© 2013 by the American Congress of Rehabilitation Medicine

A reduction in the opening diameter of the upper esophageal sphincter (UES) during the pharyngeal stage of swallowing is a physiological impairment underlying dysphagia in both neurologic and oncologic etiologies. Mechanical and neurologic events that facilitate UES dilation during the pharyngeal swallow include central nervous system—mediated relaxation, bolus pressure at the superior UES, and traction applied to UES tissue via superior and anterior movement of the hyoid and thyroid (hereafter referred to as hyolaryngeal excursion) secondary to activation in the hyolaryngeal (mylohyoid, geniohyoid, anterior digastric, and thyro-hyoid) musculature.<sup>1-5</sup> Reduced hyolaryngeal excursion is a common cause of aspiration in dysphagic patients.<sup>6,7</sup> When hyolaryngeal excursion is restricted, the diminished superior and anterior laryngeal movement fails to provide adequate traction on

0003-9993/13/36 - see front matter © 2013 by the American Congress of Rehabilitation Medicine http://dx.doi.org/10.1016/j.apmr.2013.04.013

No commercial party having a direct financial interest in the results of the research supporting this article has conferred or will confer a benefit on the authors or on any organization with which the authors are associated.

Surface electrodes and chin braces used in this study were donated to the Department of Communication Sciences & Disorders by Ampcare, LLC. The author has no financial interest or conflict of interest with this company.



INT J LANG COMMUN DISORD, MARCH–APRIL 2018, VOL. 53, NO. 2, 405–417

# **Research Report**

# Combined electrical stimulation and exercise for swallow rehabilitation post-stroke: a pilot randomized control trial

Lise Sproson<sup>†</sup> , Sue Pownall<sup>‡</sup>, Pam Enderby<sup>§</sup> and Jenny Freeman<sup>¶</sup>

†NIHR Devices for Dignity HTC, Sheffield Teaching Hospitals NHS Foundation Trust, Sheffield, UK
‡Speech & Language Therapy Department, Sheffield Teaching Hospitals NHS Foundation Trust, Sheffield, UK
§School of Health and Related Research, Sheffield University, Western Bank, Sheffield, UK
§Faculty of Medicine and Health, Leeds University, Leeds, UK

(Received September 2016; accepted October 2017)

### Abstract

*Background:* Dysphagia is common after stroke, affecting up to 50% of patients initially. It can lead to poststroke pneumonia, which causes 30% of stroke-related deaths, a longer hospital stay and poorer health outcomes. Dysphagia care post-stroke generally focuses on the management of symptoms, via modified oral intake textures and adapted posture, rather than direct physical rehabilitation of the swallowing function. Transcutaneous neuromuscular electrical stimulation (NMES) is a promising rehabilitation technology that can be used to stimulate swallowing; however, findings regarding efficacy have been conflicting.

*Aims:* This pilot randomized controlled study involving three UK sites compared the efficacy of the Ampcare Effective Swallowing Protocol (ESP), combining NMES with swallow-strengthening exercises, with usual care in order to clarify evidence on NMES in the treatment of dysphagia post-stroke. A further objective was to pilot recruitment procedures and outcome measures in order to inform the design of a full-scale trial.

*Methods & Procedures:* Thirty patients were recruited and randomized into either (1) usual speech and language therapy dysphagia care; or (2) Ampcare ESP, receiving treatment 5 days/week for 4 weeks. Outcome measures included: the Functional Oral Intake Scale (FOIS), the Rosenbek Penetration-Aspiration Scale (PAS) and patient-reported outcomes (Swallow Related Quality of Life—SWAL-QOL).

*Outcomes & Results:* Thirty patients were recruited; 15 were randomized to the Ampcare ESP intervention arm and 15 to usual care. A greater proportion (75%, or 9/12) of patients receiving Ampcare ESP improved compared with 57% (or 8/14) of the usual-care group. Patients receiving Ampcare ESP also made clinically meaningful change (a comparative benefit of 1.5 on the FOIS, and on the PAS: 1.35 for diet and 0.3 for fluids) compared with usual care. The intervention group also reported much better outcome satisfaction.

*Conclusions & Implications:* The pilot demonstrated successful recruitment, treatment safety and tolerability and clinically meaningful outcome improvements, justifying progression to a fully powered study. It also showed clinically meaningful treatment trends for the Ampcare ESP intervention.

Keywords: dysphagia, stroke, electrical stimulation, randomized controlled trial, rehabilitation.

### What this paper adds

What is already known on the subject

There is a growing movement to identify dysphagia interventions that can restore swallow function rather than simply manage symptoms. One method under evaluation is NMES; however, research to date on its efficacy has yielded conflicting results, although there is a growing consensus on its benefit as an adjunct to therapy. Therefore, we conducted a trial to work towards greater clarity to inform clinical practice.

Address correspondence to: Lise Sproson, NIHR Devices for Dignity HTC, Sheffield Teaching Hospitals NHS Foundation Trust, Sheffield S10 2JF, UK; e-mail: Lise.Sproson@nhs.net

International Journal of Language & Communication Disorders ISSN 1368-2822 print/ISSN 1460-6984 online © 2017 Royal College of Speech and Language Therapists DOI: 10.1111/1460-6984.12359

### **RESEARCH ARTICLE**

**Open Access** 



# The effect of transcutaneous neuromuscular electrical stimulation on laryngeal vestibule closure timing in swallowing

Christopher R. Watts<sup>\*</sup> and Matthew J. Dumican

### Abstract

**Background:** The purpose of this study was to investigate the effect of transcutaneous neuromuscular electrical stimulation (NMES) on the timing of laryngeal vestibule closure during the pharyngeal stage of swallowing in healthy adults. The theoretical framework proposed that NMES applied to these muscles would present a perturbation to laryngeal vestibular closure reaction time (the amount of time for the laryngeal vestibule to close once the swallowing reflex has been triggered) by providing an antagonistic force to the direction of vestibule closure.

**Methods:** Nine healthy adults (2 males, 7 females) received ten consecutive stimulations applied to the submandibular hyolaryngeal muscles while performing dry swallows. Laryngeal vestibule closure reaction time (LVCrt) and the laryngeal vestibule closure duration (LVCd) were measured from videoflouroscopic images pre-stimulation and post-stimulation.

**Results:** Results indicated a significant effect of stimulation on LVCrt but not LVCd. LVCrt was significantly reduced (timing was faster) during swallows immediately after stimulation compared to pre-stimulation.

**Conclusions:** Findings from this study support the supposition that laryngeal muscles respond to perturbations via adaptation learning, which might be used for rehabilitation of neuromuscular swallowing impairment. This pilot study supports the need for further research.

Keywords: Swallowing, Deglutition, Deglutition disorders, Neuromuscular electrical stimulation

### Background

Neuromuscular electrical stimulation (NMES) is used by allied health professionals as a treatment modality for neuromotor impairments including muscle weakness, muscle atrophy, and decreased range of motion [1-3]. In populations served by speech-language pathologists, these impairments are often associated with dysphagia caused by etiologies such as stroke and degenerative disease. The clinical research literature reporting outcomes for the effects of NMES on swallowing function is growing, as are the stimulation devices and NMES clinical protocol options available to professionals. The physiological adaptations secondary to NMES application for muscles involved during swallowing, specifically muscle

\* Correspondence: c.watts@tcu.edu



Clinical goals when utilizing NMES often include strengthening, increasing tone (to reduce atrophy), or increasing range of motion in a target muscle. In order to increase muscular strength and tone, treatments often overload (stress) a muscle by providing a resistance to contraction force. When a resistance is applied against contraction, the muscle is stressed and over time the body will respond by adapting to the stress. One way muscles adapt to the stress of an overload or resistance against contraction is to hypertrophy. Muscle hypertrophy results from the muscle cells gaining volume.

© The Author(s). 2018 **Open Access** This article is distributed under the terms of the Creative Commons Attribution 4.0 International License (http://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made. The Creative Commons Public Domain Dedication waiver (http://creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated.

Davies School Communication Sciences & Disorders, Texas Christian University, TCU Box 297450, Fort Worth, TX 76129, USA



Article

# Neuromuscular Electrical Stimulation Plus **Rehabilitative Exercise as a Treatment for Dysphagia** in Stroke and Non-Stroke Patients in an NHS Setting: **Feasibility and Outcomes**

### Nicola Martindale <sup>1,\*</sup>, John Stephenson <sup>2</sup>, and Sue Pownall <sup>1,\*</sup>

- 1 Sheffield Teaching Hospitals NHS Foundation Trust; Royal Hallamshire Hospital, Glossop Road, Sheffield S10 2JF, UK
- 2 University of Huddersfield, Queensgate, Huddersfield HD1 3DH, UK; J.Stephenson@hud.ac.uk
- Correspondence: Nicola.martindale@nhs.net (N.M.); sue.pownall2@nhs.net (S.P.)

Received: 13 August 2019; Accepted: 17 September 2019; Published: 24 September 2019



Abstract: Dysphagia is a debilitating condition with significant consequences in terms of physical and mental health. This study demonstrates that it is feasible to provide an intensive therapy program combining neuromuscular electrical stimulation (NMES) with exercise against resistance in the treatment of dysphagia in a public healthcare setting. Thirty-one patients (17 stroke, 14 non-stroke) who experienced dysphagia with reduced laryngeal elevation completed the therapy program. After checking the data sets for comparability, it was deemed appropriate for the outcome data from these patients to be combined with that of 12 stroke patients previously reported to enable statistical analysis on a larger data set (n = 43). A repeated-measures ANOVA revealed a statistically significant increase in amount and variety of food a patient was able to take orally (FOIS) following completion of treatment (p < 0.001). There was no significant between-subject effect of stroke status (p = 0.43), or interaction between treatment and stroke status (p = 0.68). There was a significant improvement in secondary outcome measures of swallow safety with fluids (PAS) (p < 0.001) and swallow-related quality of life (Swal-Qol (p < 0.001). These findings indicate that the therapy program may be associated with reduced impairment in a subset of patients with dysphagia resulting from stroke and non-stroke atiologies, and the data will inform the design of future research to address unanswered questions.

Keywords: dysphagia; stroke; rehabilitation; neuromuscular electrical stimulation; NMES; Ampcare ESP

### 1. Introduction

The use of neuromuscular electrical stimulation (NMES) in the treatment of dysphagia has received considerable research interest (see the National Institute for Clinical Excellence (NICE) guidelines, 2018, for an overview [1]). Dysphagia, or difficulty moving material from the mouth to the stomach, is a debilitating condition with far-reaching consequences in terms of a patient's physical health [2,3], mental wellbeing [4], level of dependency and economic cost to the state [2,5]. Dysphagia can have many causes, including stroke, neurodegenerative conditions, traumatic brain injury and head and neck cancers. Severity of dysphagia also varies widely, from some people having a mild difficulty that may result in the avoidance of certain foods or eating situations; to those with severe cases, rendering a patient unable to safely eat or drink orally at all, thus requiring consideration of clinically assisted artificial nutrition and hydration. Standard management of dysphagia in the UK may involve exercise regimens, but is often restricted to symptom management through modification of oral intake and Geriatrics 2019, 4, 53; doi:10.3390/geriatrics4040053 www.mdpi.com/journal/geriatrics



### Ampcare ESP in Treatment of Post-Stroke Dysphagia Further Evidence for Positive Outcomes

IO<sup>™</sup> ESSD <u>CONGRESS</u> Dysphagia Management across Lifespa New State of the Art post Covid-19? 8-10 OCT 2020

Neuromuscular electrical stimulation plus intensive exercise against resistance (Ampcare ESP) in the treatment of post-stroke dysphagia: further evidence for positive outcomes

Dysphagia in respiratory diseases / stroke and brain damage **PB08** 

### Introduction

Dysphagia is a common consequence of stroke. It can result in malnutrition and dehydration, aspiration and pneumonia, and reduced quality of life. There are also significant financial costs incurred in the management of dysphagia and associated consequences post stroke, further motivating the search for effective treatments. There is evidence which indicates that a therapy program combining neuromuscular electrical stimulation (NMES) to the suprahyoid muscles with the search content contents on the fractive Effective Surallowing Pertocol, can be intensive exercise against resistance (the Ampcare Effective Swallowing Protocol) can be effective in the treatment of post stroke dysphagia<sup>1</sup>. More data is required to further evaluate the effectiveness of the therapy program for this patient group. In accordance with NICE guidelines<sup>2</sup> outcome measures are collected in the UK when Ampcare ESP is used in the treatment of dysphagia. The outcomes from 26 stroke patients who were offered the therapy program across SUIS to use the factore data is insurand.

### Objective

Collate outcome data from across the UK to determine what proportion of patients experience a functional gain in their swallow in terms of improved quantity and or variety of oral intake following completion of Ampcare ESP.

### Methods

26 stroke patients from 5 NHS trusts in the UK, identified by Speech and Language Therapists as having pharyngeal dysphagia including reduced laryngeal elevation (either during bedside swallow assessment or videofluoroscopy) were offered up to 22 thirty minute sessions of Ampcare ESP

Ampcare ESP = NMES to suprahyoid muscles + exercise against resistance

Functional oral intake scale (FOIS)<sup>3</sup> was recorded before and after completion of the therapy program.

Of the 26 patients offered Ampcare ESP, 3 died and one was too unwell to complete the treatment.



22 patients completed the therapy program (13 male, age range 31-91, mean age 69 years) 19 (86%) were able to safely increase the amount/ variety of food they were eating 5 SIO 4 = total oral diet 2 = minimal oral intake Ш ш 1 = NBM Patients FOIS before therapy

Results

10/11 patients who were NBM at the start of therapy were able to introduce some oral intake 9/20 patients who were tube dependent were able to meet their needs orally after therapy

### Conclusions

- Speech and language therapists in NHS services across the UK are using Ampcare ESP in the treatment of post stroke dysphagia
- Outcome data continue to demonstrate that many of these patients are able to eat a more varied diet following completion of the therapy program.
- A large proportion of these patients were able to meet their needs orally after completing the therapy program.
- Further research is required to evaluate the health economics of the therapy program.

References: 1. Martindale et al., 2019 Geriatrics, 4, 53; 2. NICE https://www.nice.org.uk/guidance/ipg634; 3. Crary et al., 2015 Arch. Phys. Med. Rehabil. 2005, 86, 1516-1520

### Is There a Role for Ampcare ESP in the Treatment of Post Covid-19 Dysphagia?

OTH ESSD CONGRESS Dysphagia Managem New State-of-the-Art p 8-10 OCT 2020

Is there a role for neuromuscular electrical stimulation therapy programs (e.g. the Ampcare Effective Swallowing Protocol) in the treatment of post covid-19 dysphagia? Nicola Martindale, Amanda Forrester\*, Jude Sadiq and Sue Pownall Sheffield Teaching Hospitals NHS Foundation Trust UK (\*Amanda.forrester1@nhs.net)

Dysphagia during/ post covid-19 **PD08** 

### Introduction

A cohort of four patients with Covid-19 requiring a lengthy critical care stay including intubation, ventilation and tracheostomy presented with severe pharyngeal dysphagia. Following decannulation, videofluoroscopy (VF) analysis showed reduced pressure generation, reduced laryngeal elevation & closure, reduced base of tongue to posterior wall contact & absent epiglottic inversion leading to reduced pharyngeal clearance and penetration' aspiration. All patients were unsafe to commence oral intake and were artificially fed. Their dysphagia did not improve with standard swallow therapy. The Ampcare Effective Swallow Protocol (ESP) combines neuromuscular electrical stimulation (NMES) with swallow exercises completed against resistance. Evidence suggests that this may be beneficial in treating some patients with pharyngeal dysphagia'. A trial of Ampcare ESP was introduced with these patients to determine whether they were able to tolerate and actively participate with the therapy whilst recovering from Covid-19 and whether they made functional gains to their swallow.

### Objectives

Are post-critical care Covid 19 patients able to tolerate and engage with Ampcare ESP? Do patients see an improvement in their functional swallow after engaging with Ampcare ESP?

### Methods

The following four covid-19 survivors with severe pharyngeal dysphagia were offered Ampcare ESP						
	Age	Gender	Intubation (days)	Tracheostomy (days)	Dysphagla onset – Ampcare ESP (days)	
А	73	м	11	41	65	
В	56	М	10 + 6 (intubated twice)	25	47	
С	66	М	9	37	77	
D	67	М	32	20	57	

Each patient was offered 20 thirty minute sessions of Ampcare ESP

Ampcare ESP = NMES to suprahyoid muscles + exercise against resistance

Swallow was assessed before and after treatment at bedside and with VF (\*not all VF had been completed at time of writing)

Outcome measures were the penetration and aspiration scale (PAS)<sup>2</sup> and Waxman dysphagia severity score<sup>3</sup> as measured at VF and level of oral intake.

Results							
	Sessions completed	Average no. of therapy sessions per week	PAS pre PAS post	Waxman pre Waxman post	Oral intake pre Oral intake post		
A	20/20	3.3	3-5 5	6 5	NBM Trials of level 5 diet		
В	17/20	2.4	3-5 3-5	6 5	NBM Trials of level 5 diet		
С	20/20	2.8	7 7	5 5	NBM Trials of level 5 diet		
D	15/20	3	5 *	6 *	NBM Trials of level 5 diet		

All four patients tolerated the therapy well & were able to reliably generate swallows in synch with stimulation Some improvement in swallow was observed following treatment in terms of reduction in Waxman dysphagia severity score and safe introduction of trials of oral intake (diet) post therapy

Patients continued to experience ongoing moderately-severe dysphagia despite full engagement with Ampcare ESP

All patients continued to penetrate/ aspirate fluids. Patient A had a further VF 6 weeks after completion of Ampcare ESP, their swallow had improved further, allowing safe introduction of small signs of level 0 fluids and level 6 diet

### Conclusions

- This cohort of post critical care Covid-19 patients tolerated & were motivated to engage with Ampcare ESP Their dysphagia was severe in presentation but improved sufficiently post therapy to allow introduction of trials of level 5 diet. Unfortunately it was not possible for patients to safely commence fluids immediately after completing therapy.
- The dysphagia presentation was not typical of the disuse atrophy & deconditioning seen in standard post-critical The dysphagia presentation was not typical of the disuse atrophy & deconditioning seen in standard post-critical
- The opporting presentation was not plocate on the based exploring determination of the plocate of plocate of plocate of the pl

References: 1. Martindale et al., 2019 Geriatrics, 4, 53; 2. Rosenbek et al., 1996, Dysphagia, 11, 93-98. 3. Waxman et al., 1990, Nutr Clin Pract, 5:196-199
#### ORIGINAL ARTICLE

## Immediate Effects of Electrical Stimulation on **Oropharyngeal Structure and Laryngeal Vestibular Closure: A Pilot Study in Healthy Subjects**

Michihiro Ogura, SLHT <sup>a</sup> Shuji Matsumoto, MD, PhD <sup>b,c</sup> Rintaro Ohama, MD <sup>c</sup> Yumi Ohama, MD <sup>c,d</sup> Haruka Arima, SLHT <sup>a</sup> Keita Takenaka, SLHT <sup>e</sup> Keiichi Toyama, SLHT <sup>f</sup> Toshiyuki Ikegami, SLHT <sup>g</sup> and Megumi Shimodozono, MD, PhD <sup>c</sup>

Objectives: This study examined the immediate effects of neuromuscular electrical stimulation (NMES) on the dynamics of oropharyngeal structure and laryngeal vestibular closure (LVC) in healthy subjects. Methods: Ten healthy male volunteers participated in this controlled, beforeand-after, videofluoroscopic swallowing pilot study. The study was conducted in four phases (each performed twice): (1) saliva swallow (SS) before evaluation (BEFORE), (2) NMES while at rest with no SS (NMES AT REST), (3) SS during NMES (DURING NMES), and (4) SS to examine the aftereffects of NMES (AFTER). We measured distances that oropharyngeal structures moved in the NMES AT REST phase, and we analyzed the kinematics of saliva swallowing primarily in the BEFORE and AFTER phases. **Results:** Four changes in the morphology of the oropharyngeal structure caused by NMES AT REST were statistically significant: anterior-upward displacement of the hyoid bone and larynx, stretch of the laryngeal vestibule, and posterior ridge of the tongue root. Regarding the kinematics measured during SS, although there was no significant change in LVC reaction times, LVC duration in the AFTER phase was significantly longer than BEFORE. Regarding maximal displacement of the hyoid bone, there was significantly greater movement AFTER than BEFORE. As additional exploratory outcomes, the velocity of hyoid bone movement was significantly slower, and the hyoid-to-larynx approximation was significantly smaller, DURING NMES than AFTER. Conclusions: Longer duration of LVC might be caused by adaptive learning with NMES-induced structural changes in the oropharynx. Further clinical studies are warranted to determine whether this approach improves dysphagia, which impairs LVC.

**Key Words:** Ampcare ESP<sup>TM</sup>; hyoid laryngeal proximity; laryngeal vestibular closure; neuromuscular electrical stimulation; swallowing

#### INTRODUCTION

nism preventing penetration or aspiration during swallowing. Given that impaired LVC causes unsafe swallowing,

Laryngeal vestibule closure (LVC) is the primary mecha-

- <sup>a</sup> Department of Rehabilitation, Kagoshima University Hospital, Kagoshima, Japan
- <sup>b</sup> Department of Rehabilitation and Physical Medicine, Mito Clinical Education and Training Center, Tsukuba University Hospital, Mito-kyodo Hospital, Mito, Japan
- ° Department of Rehabilitation and Physical Medicine, Kagoshima University Graduate School of Medical and Dental Sciences, Kagoshima, Japan
- <sup>d</sup> Neurosurgery Center, Ichikikushikino Medical Association, Kagoshima, Japan
- e Department of Rehabilitation, Tarumizu Municipal Medical Center, Tarumizu Chuo Hospital, Kagoshima, Japan
- <sup>f</sup> Department of Speech Language-Hearing Therapy, Faculty of Allied Health Sciences, Yamato University, Osaka, Japan
- <sup>g</sup> Department of Rehabilitation, Kikuno Hospital, Kagoshima, Japan

Copyright © 2022 The Japanese Association of Rehabilitation Medicine



This is an open-access article distributed under the terms of the Creative Commons Attribution Non-Commercial No Derivatives (CC BY-NC-ND) 4.0 License. http://creativecommons.org/licenses/by-nc-nd/4.0/

Received: November 8, 2021, Accepted: June 7, 2022, Published online: July 9, 2022

Correspondence: Rintaro Ohama, MD, 8-35-1 Sakuragaoka, Kagoshima City, Kagoshima 890-8520, Japan, E-mail: rintaro@m2.kufm. kagoshima-u.ac.jp

Swallowing Out	Comes in Clinically Based Pat Stimulation ( Kaitlynn Harper Matthew Dumican, Ph.	ients Using Neuromuscular VMES) <sup>B.S.</sup> O. CCC-SLP	Electrical Mestern Michigan
<ul> <li>Purpose</li> <li>Purpose</li> <li>Neuromuscular electrical stimulation (NMES) has received mixed findings in the literature regarding its effectiveness in treating dysphagia of neurogenic etiology (Speyer et al., 2022). While functional outcomes are often considered, the underlying mechanisms of improvement are often unspecific, examining outcomes such as water swallow tests, "swallowing ability", or only penetration or aspiration (Alamer et al., 2020, Bath et al., 2018).</li> <li>The aims of this project were to examine the effects of NMES on mechanistic changes to swallow tests, "swallowing ability", or only penetration or aspiration in:</li> <li>The aims of this project were to examine the effects of NMES on mechanistic changes to swallow tests, "swallow ability and movement distance 2. Laryngeal kinematic timings 3. Swallow safety outcomes</li> <li>Methods</li> <li>Comparison of baseline and final VFSS of 16 (13 mole, 3 female) of clinical patients with post-stroke after a range of 73-5 days a week and a minimum of 30 total days of treatment to obtain:</li> <li>Hyoid Kinematics (anterior, superior, excursion)</li> <li>Time-To-Laryngeal Vestibule Closure (LVC) &amp; LVC Duration</li> <li>Penetration Stander (VFSS review. All patients included in this study received NMES as a treatment protocol for dysphagia relabilitation as their primary treatment, in various clinical settings. NMES was implemented with suprativoid placement, pulse duration of 50 µs, and variable internetive in macions clinical settings.</li> </ul>	Pre-Treatment Typical and Atypical Swallow Safety (Figure 1)	Post-Treatment Typical and Atypical Swallow Safety (Figure 2)	<b>Conclusions &amp;</b> <b>Clinical</b> <b>Relevance</b> The results from this project support that NMES may be a viable treatment option for patients experiencing post- stroke dysphagia, in agreement with findings from Speyer et al. (2022), and particularly in dysphagia characterized by reduced laryngeal elevation and events of penetration and/or aspiration (Martindale et al., 2019). In a small, clinically based cohort of patients with post- stroke dysphagia, NMES with suprahyoid placement and parameters based on muscle and swallow physiology significantly improved hyoid excursion, decreased inte-lo- LVC, and decreased occurrences of penetration and/or aspiration. More rigorously controlled pre- post treatment studies are required to support NMES treatment efficacy.
and progression. A one-way repeated measures ANOVA (RM- MANOVA) was performed to examine within- subjects changes of hyoid kinematics and laryngeal vestibule timings pre- and post-treatment. Crosstabulations were performed to examine pre- and post-treatment associations with penetration and/or aspiration.	Time-to-LVC Post-Treatment Average LVC Post-Treatment	of treatment on hyoid excursion (F= 9.32, $p=$ 0.008) and time-to-LVC (F= 54.74, $p<$ 0.001) (figure 3). Effect size calculations of partial eta-squared ( $\eta^2_p$ ) indicated that both hyoid excursion ( $\eta^2_p$ = .37) and time-to-LVC ( $\eta^2_p$ = .77) demonstrated large effects.	References

## **Strategies to Increase Swallowing Activities**

#### Good oral care should be maintained before therapy

- 1. Cold, wet spoon: present spoon with just cold condensation to initiate the autonomic response
- 2. Lollipops (sugar-free if diabetic) or ice lollies, laryngeal mirrors are good options to stimulate sensory
- receptors on the mucosa of the mouth, tongue, and pharyngeal regions to trigger a swallow
- 3. Alterations in bolus volume, taste, temperature
- 4. Taste Tests, using food coloring or hidden juice cups; appropriate consistencies; & diabetic considerations
- 5. Smell or taste of lemon, cranberry, apple, or grape juices to increase saliva
- 6. Coconut oil for dry tongue/xerostomia to assist and lubricate the surface of the tongue
- 7. Consider a water protocol, with appropriate oral care

#### **Pharyngeal Exercises and Strategies**

#### Valleculae/Tongue Base

- 1. Tongue slides
- 2. Yawn
- 3. Masako (protrude tongue between teeth hold swallow hard)
- 3. Pretend to gargle
- 4. Pull/hump tongue base posteriorly
- 5. Head tilt or turn as appropriate
- 6. Resistive Chin-to-Chest for 5 seconds with RPD and effortful swallow
- /k/ and /g/ words/sounds (say "kuh" or "guh")
- 8. Mouthguard pull

# PES/UES/Cricopharyngeal sphincter if related to reduced laryngeal elevation

- 1. Effortful/dry/hard swallow
- 2. Alternate solids and liquids if appropriate
- 3. Head turn, lift, or tilt if appropriate
- 4. Chin thrust (not chin tuck!)
- 5. Hyolaryngeal exercises /Mendelsohn Maneuver

#### **Pharyngeal Pyriform**

- 1. Head turn to impaired side if appropriate
- 2. Tilt head to stronger side if appropriate
- 3. Alternate solid and liquids if appropriate
- 4. Push against mandible while swallowing
- 5. Successive swallow

#### Pharyngeal Posterior wall

- 1. Repeat effortful swallows and increase rate of swallow
- 2. Thermal/Tactile/Chemical stimulation to elicit effortful swallows

#### Laryngeal or Hyolaryngeal exercises

- 1. Resistive Chin-to-Chest for 5 seconds with RPD & effortful swallow
- 2. Resistive Chin-to-Chest with RPD & Mendelsohn Maneuver for 5 seconds
- Resistive Chin-to-Chest for 5 seconds with RPD & jaw open/close then effortful swallow
- Shaker exercise: Lie flat raise head look at toes/feet - hold for 1 minute then relax for 1 minute, repeat 2 more times, followed by 30 consecutive head raises while keeping shoulders resting flat.
- 5. EMST
- 6. Falsetto exercises
- 7. Pitch glides gliding up in pitch to reach a high, squeaky voice, then holding it for a minimum of 5 seconds with effort.
- 8. Adduction exercises while lifting up or pushing down on chair
- 9. Supraglottic swallow
- 10. Valsalva maneuver (breath hold) contraindicated with cardiac precautions

#### Home-exercise program for Head & Neck Cancer patients – 5 sets of 5 reps, 4 times per day

- 1. Tongue press
- 2. Effortful swallows
- 3. Falsettos (effortful "eeeee" pronunciation)
- 4. Chin-to-Chest with jaw open against resistance (RPD)
- 5. While sitting, tilt head back with jaw open attempt to close jaw and protrude it toward the ceiling until a stretch is felt on the front of the neck. Hold this for 30 seconds for 5 reps.

## **Oral Motor Function**

#### 5 sets of 5 reps, 4 times per day

Lips – reduced lip sensation, strength, and ROM may result in drooling, lip biting and pocketing, as well as anterior loss.

- 1. Assess labial musculature by pressing lips tightly around a spoon, mouthguard, or tongue depressor as the SLP tries to remove it.
- 2. Straw sucking against resistance. Attempt holding stickers, coins or other weighted objects at end of straw.
- 3. Have patient close eyes. Lightly touch each of the 4 quadrants of the upper and lower lips with a tongue depressor or cotton swab. This will assess lip sensitivity.
- 4. Smile/Pucker. Have them say "ee", Have them say "oo" / Alternate between "ee" and "oo".
- 5. Say the syllable "puh" as quickly as possible to determine ability to obtain closure of lips.
- 6. Open mouth but try to hold a pucker.
- 7. Blowing Bubbles/Blowing out a Candle; Horn/Harmonica/Party Favor; Straw– controlled pursing of the lips.
- 8. Close lips around the tongue blade/Press lips together "MMMMMMM".

**Cheeks** –reduced cheek sensation, strength, and range of motion may result in pocketing of food or biting cheek. Before beginning exercises, observe patient's facial asymmetry. Check buccal musculature by:

- 1. Have them puff their cheeks up with air. Be sure to check velum function before this exercise because reduced velopharyngeal seal can result in decreased oral pressure and may contribute to the patient's inability to puff cheeks.
- 2. Assess the sensitivity of cheeks by touching random locations with a cotton swab.
- 3. Pretend to "swish mouthwash around in the mouth".
- 4. Push cheeks out with spoon while the patient attempts to suck cheeks in.

**Tongue** – reduced tongue sensation, strength, and ROM can result in the inability to manipulate bolus in the oral cavity, it can also result in residue on hard palate and pocketing.

- 1. Assess sensitivity by touching various areas of the tongue.
- 2. Stick out tongue check for deviation, fatigue, and fasciculations
- 3. Move tongue from left to right check for ROM, speed, and coordination.
- 4. Elevate tip of tongue maintain for 5 seconds. Can also elevate tip to alveolar ridge and move it anteriorly to posteriorly along palate (tongue slide).
- 5. Using Peanut Butter on lips to increase side-to-side ROM.
- 6. Stick out tongue and point down to chin maintain for 5 seconds, left to right cheek, maintain 5 seconds.
- 7. Tongue retraction humping it posteriorly and hold for 5 seconds, /g/ and /k/ words/sounds/soundless "guh.".
- 8. Lingual resistance Extend tongue forward and push against top row of teeth/palate, tongue depressor, lollipop, or finger. Hold for 5 seconds.
- 9. Diadochokinetic rate: "Buttercup".

## **Soft Palate** – reduced sensation, strength, and ROM can result in reduced gag reflex, hypernasality, nasal reflux, and premature spillage

- 1. Phonate "ah" in isolation and repetitive pattern observe elevation.
- 2. Assess for hypernasality.
- 3. Assess gag reflex contraindicated with cardiac precautions.

#### Velar Exercises

- 1. Thermal Tactile Stimulation to help elicit effortful swallowing exercise.
- 2. Suck/swallow; Straw attempts with Thick Liquid to increase effort; Can also suck through pinched straw.
- 3. Effortful "guh".
- 4. High pitch "eeee" for 5 seconds.
- 5. Snore for 5 seconds.

Jaw strength and ROM – Necessary for rotary movements involved in mastication. This can be used with patients exhibiting reduced oral ROM or masseter weakness.

Mandible - mastication/ side-to-side and rotary action is required for cohesive boluses

- 1. Open mouth wide, hold for 5 seconds.
- 2. Open and close 5 times quickly.
- 3. Mouthguard pulls.
- 4. Jaw grading (bite and hold for 5 seconds) using chewy tubes.
- 5. Move jaw from left to right and in a circular motion as though rotary chewing.

#### www.ampcarellc.com



900 42<sup>nd</sup> Street South PO Box 6757 Fargo, ND 58108-6757 PDAC Medicare Pricing, Data Analysis and Coding

February 24, 2009

#### LERKINS BILLING & CONSULTING PROFESSIONALS INC ATTN: CHRYSTAL ZELLER 407 E MARKET ST STE 105B CRAWFORDSVILLE IN 47933

Re: Restorative Medical's RPD (Restorative Posture Device) Model #50707

Xref #: 7493545

Dear Ms. Zeller:

This letter is in response to your recent inquiry for coding verification of the above listed product(s) manufactured and/or distributed by your company.

It is our determination that the Medicare HCPCS code(s) to bill the four DME Medicare Administrative Contractors (DME MACs) is/are:

L0140 - CERVICAL, SEMI-RIGID, ADJUSTABLE (PLASTIC COLLAR)

This HCPCS coding decision applies to the submitted product(s) as presented to and reviewed by the Pricing, Data Analysis and Coding (PDAC). Any modifications to the product(s) could change the HCPCS code and would need to be reviewed for coding verification. The assignment of a HCPCS code to the product(s) should in no way be construed as an approval or endorsement of the product(s) by PDAC or Medicare, nor does it imply or guarantee claim reimbursement or coverage. For questions regarding claim coverage or reimbursement, please contact your DME MAC jurisdiction.

Should you disagree with this coding decision, a reconsideration can be requested. The reconsideration form is located on the PDAC website (<u>www.dmepdac.com</u>) under the HCPCS Review tab. The PDAC will reconsider the request if made within 45 days of the date of this letter and additional supporting documentation is provided. If a request for reconsideration is made after 45 days, the request is treated as a new Coding Verification Review and a complete application must be submitted, along with the additional documentation supporting the request.



www.ampcarellc.com

#### HCPCS Code: L0140 Short Description: Cervical semi-rigid adjustab Long Description: CERVICAL, SEMI-RIGID, ADJUSTABLE (PLASTIC COLLAR) https://www4.palmettogba.com/pdac\_dmecs/ hcpcsdetailsFeeScheduleSearch.do?hcpcsCode=L0140

Beneficiary State	Non-Rural	Effective		Beneficiary State of	Non-Rural	Effective	
of Residence	Fee	From	Effective To	Residence	Fee	From	Effective To
AL	\$82.58	1/1/2024	3/31/2024	NV	\$70.43	1/1/2024	3/31/2024
AK	\$61.46	1/1/2024	3/31/2024	NH	\$91.05	1/1/2024	3/31/2024
AZ	\$70.43	1/1/2024	3/31/2024	NJ	\$70.43	1/1/2024	3/31/2024
AR	\$71.86	1/1/2024	3/31/2024	NM	\$71.86	1/1/2024	3/31/2024
CA	\$70.43	1/1/2024	3/31/2024	NY	\$70.43	1/1/2024	3/31/2024
CO	\$79.52	1/1/2024	3/31/2024	NC	\$82.58	1/1/2024	3/31/2024
СТ	\$91.05	1/1/2024	3/31/2024	ND	\$79.52	1/1/2024	3/31/2024
DE	\$71.91	1/1/2024	3/31/2024	ОН	\$70.43	1/1/2024	3/31/2024
DC	\$71.91	1/1/2024	3/31/2024	ОК	\$71.86	1/1/2024	3/31/2024
FL	\$82.58	1/1/2024	3/31/2024	OR	\$93.91	1/1/2024	3/31/2024
GA	\$82.58	1/1/2024	3/31/2024	PA	\$71.91	1/1/2024	3/31/2024
HI	\$65.74	1/1/2024	3/31/2024	PR	\$136.68	1/1/2024	3/31/2024
ID	\$93.91	1/1/2024	3/31/2024	RI	\$91.05	1/1/2024	3/31/2024
IL	\$70.43	1/1/2024	3/31/2024	SC	\$82.58	1/1/2024	3/31/2024
IN	\$70.43	1/1/2024	3/31/2024	SD	\$79.52	1/1/2024	3/31/2024
IA	\$93.91	1/1/2024	3/31/2024	TN	\$82.58	1/1/2024	3/31/2024
KS	\$93.91	1/1/2024	3/31/2024	ТХ	\$71.86	1/1/2024	3/31/2024
KY	\$82.58	1/1/2024	3/31/2024	UT	\$79.52	1/1/2024	3/31/2024
LA	\$71.86	1/1/2024	3/31/2024	VT	\$91.05	1/1/2024	3/31/2024
ME	\$91.05	1/1/2024	3/31/2024	VI	\$70.43	1/1/2024	3/31/2024
MD	\$71.91	1/1/2024	3/31/2024	VA	\$71.91	1/1/2024	3/31/2024
MA	\$91.05	1/1/2024	3/31/2024	WA	\$93.91	1/1/2024	3/31/2024
MI	\$70.43	1/1/2024	3/31/2024	WV	\$71.91	1/1/2024	3/31/2024
MN	\$70.43	1/1/2024	3/31/2024	WI	\$70.43	1/1/2024	3/31/2024
MS	\$82.58	1/1/2024	3/31/2024	WY	\$79.52	1/1/2024	3/31/2024
MO	\$93.91	1/1/2024	3/31/2024				
MT	\$79.52	1/1/2024	3/31/2024				
NE	\$93.91	1/1/2024	3/31/2024				



Addresses anatomical alignment, posture, and strength that allows for optimal laryngeal and pharyngeal function during swallowing.

#### Indications:

The Restorative Posture Device (RPD) is indicated for patients that present with poor positioning including forward head and neck posture resulting in decreased range of motion (ROM) and increased dysphagia risk as a result of posture as determined on an instrumental exam. This patient may be **flaccid** (can be passively ranged to an upright position but cannot sustain the position), may present with **shortened adaptive tissue** (unable to be passively or actively ranged to upright position), or may present with **neurological tone** with or without shortened adaptive tissue.

#### **Intended Uses:**

Patients with impaired ability to swallow and lost cervical range of motion - with or without neurological tone – can improve function and gradually increase flexibility by reeducating muscles in the submandibular region to promote the ideal head positioning for anatomical alignment that allows for optimal hyolaryngeal excursion. Patients may also benefit from:

- 1. Enhanced swallowing treatment when used in conjunction with Ampcare's Effective Swallowing Protocol (ESP<sup>TM</sup>); assists with proper position of head and neck, maintains placement and improves electrical conduction of E Series Electrodes, and provides a resistive exercise protocol.
- 2. Improved respiratory function, coordination of respiration and swallowing and increased O<sub>2</sub> and CO<sub>2</sub> exchange.
- 3. Improved pharyngeal space for adequate laryngeal movement and passage of the bolus through the pharyngeal phase.
- 4. Increased ability to orally accept food taken by mouth (PO) and improve oral ROM, manipulation of the bolus for improved safety of the swallow.
- 5. Increased healing of any skin breakdown related to pressure, skin-on-skin contact, and lack of air to involved tissue.
- 6. Decreased neck pain and discomfort; be able to improve the ability to take in nutrition and hydration.
- 7. Improved ability to receive the stimuli of their environment which could increase their quality of life as they are able to take part in activities which may also improve their degree of orientation/cognition and PO intake.

#### **Application:**

- 1. Apply and remove per physician's order and patient's overall plan of care.
- 2. Explain to patient in understandable terms what you are about to do.
- 3. Determine the degree of comfortable range of motion of head/neck. Try to keep the cervical spine in good alignment (erect) for more appropriate vallecular space with a mild dorsal glide (chin retraction) to allow for better epiglottic closure and airway protection.
- 4. If needed this device can be remolded. Refer to the Instructions on Heat Molding section.
- 5. Adjust the blue strap around patient's neck and secure with hook-and-loop fasteners to form the most comfortable and therapeutic angle. Assure moldable chin cup is comfortable and no plastic is touching the patient's throat.
- 6. Begin wearing schedule with short duration according to the patient's tolerance levels and their individual conditions. Gradually buildup 30 minutes per day to typically no more than 6 hours per 8-hour shift (6 hours on 2 hours off). Attempt to build muscle strength and prevent weakening. Consider initiating use during ESP dysphagia treatment and/or liquid PO trials.
- 7. Check all areas of skin where the device touched each time it is removed. If any red areas appear that do not resolve within 30 minutes, remove the device and inform the caregiver or call **Ampcare at 682-561-2444**. Shortening the wearing period or adding more padding may resolve the issue.
- 8. Chart all pertinent information.

#### Laundering:

Cleaning instructions: Remove padding. Hand wash, lay flat and blot with a towel, then air dry. Wipe plastic clean with mild cleanser, protecting the hook-and-loop fasteners. Rinse with water, and allow to dry thoroughly while protecting the hook-and-loop fasteners. Replace padding, making sure all plastic edges where patient's skin might touch is covered.

## **Restorative Posture Device** SINGLE PATIENT USE ONLY

## **Instructions on Heat Molding**



- 1. Before heat-molding the brace, place the chin plate under the patient's chin to see where modifications are needed. Patient's chin should remain slightly behind the end of the chin plate.
- 2. Before applying heat with a heat gun, remove chin and chest padding, strap, and hook-and-loop fastener on S-curve.
- 3. Spot mold the device only at the middle or bottom of the S-curve. Only mold to 10-15 degrees of stretch.
- 4. Move the heat gun in a small circular motion approximately 3" away from the plastic, for approximately 30 seconds on each side, only to the area to be remolded. Heat the middle of the S-curve to adjust the height and the bottom of the S-curve to adjust the depth of the chest plate.
- 5. Replace hook-and-loop fastener, padding, and strap to appropriate areas after plastic cools.
- 6. Place the brace on the patient to ensure that desired modifications are complete.

## \*\* Not for prolonged use \*\*

You may change the curve angle, but keep the S-curve to facilitate a cushioning effect.
Spot mold only at the middle of the S-curve to adjust the height or the bottom of the S-curve to adjust the chest plate.
Avoid heating the top of the S-curve as the angle is set to the appropriate amount of resistance.

**Warranty**: Ampcare agrees that if any defect in materials or workmanship appears in this product within six months after the original date of consumer purchase, it will replace the product in question free of charge. This applies only if the product has been used as directed and has not been damaged through misuse, accident, or neglect.

If a defect appears, please check that the product is being used in accordance with instructions. If so, return it with this warranty and proof of purchase to your nearest Ampcare Dealer.

# If you are not a certified user of ESP, contact Ampcare to receive video instructions on how to perform ESP exercises.

U.S. Patent: 7,885,713

Ampcare, LLC Suite 111 1120 South Freeway Fort Worth, Texas 76104 USA Tel: (682) 561-2444 Manufactured for Ampcare by: Restorative Medical 332 East Broadway Brandenburg, Kentucky 40108 USA KSA Authorized Representative: Arabic Blue Diamond Medical Establishment Al Imam Al Shafia, Al faisaliya Jeddah 23447 Saudi Arabia



RPD-121316 Rev.C

www.ampcarellc.com

- Agency for Health Care Policy and Research. (1999). Diagnosis and treatment of swallowing disorders (dysphagia) in acute-care stroke patients. Evidence Report/Technology Assessment No. 8 (AHCPR Pub. No. 99-E024). Rockville, MD: Author.
- Ali GN, Laundl TM, Wallace KL, deCarle DJ & Cook IJS. (1996). Influence of cold stimulation on the normal pharyngeal swallow response. *Dysphagia*, 11:2-8.
- American Speech-Language-Hearing Association..(c1997-2007). Dysphagia Treatment with E-Stim: Coding Considerations.Retrieved from https://www.asha.org/Practice/reimbursement/medicare/e-stim\_coding/
- American Speech-Language-Hearing Association. (2016). Code of ethics [Ethics]. Available from www.asha.org/policy/.
- ASHA Communication Facts: Special Populations: Dysphagia- 2008 Edition. American Speech-Language-Hearing Association.
- Azola AM, Greene LR, Taylor-Kamara I, Macrae P, Anderson C & Humbert IA. (2015). The relationship between submental surface electromyography and hyo-laryngeal kinematic measures of Mendelsohn Maneuver duration. *J Speech Lang Hear Res.* 58(6):1627–36.
- Azola AM, Sunday KL, & Humbert IH (2016). Kinematic Visual Biofeedback Improves Accuracy of Learning a Swallowing Maneuver and Accuracy of Clinician Cues During Training. *Dysphagia* 32:115–122. https://doi.org/10.1007/s00455-016-9749-z
- **Baker L & Parker K.** (1986). Neuromuscular Electrical Stimulation of the Muscles Surrounding the Shoulder. *Physical Therapy*, 66(12):1930-37.
- **Baker L, McNeal D, Benton L, Bowman B & Waters R.** (1993). NeuroMuscular Electrical Stimulation, A Practical Guide, Los Amigos Research & Education Institute, 3<sup>rd</sup> Edition.
- Belafky PC, Mouadeb DA, Rees CJ, Pryor JC, Postma GN, Alen J, Leonard RJ. (2008). Validity and reliability of the Eating Assessment Tool (EAT-10). Ann Otol Rhinol Laryngol, 117(12):919-924.
- Brasseur JG & Dodds WJ. (1991). Interpretation of intraluminal manometric measurements in terms of swallowing mechanics. *Dysphagia*, 6:100-119.
- **Bickel CS, Gregory CM & Dean JC.** (2011). Motor unit recruitment during neuromuscular electrical stimulation: a critical appraisal. Eur J Appl Physiol. Oct; 111(10): 2399-407.
- **Brown M.** (2012). The physiology of age-related and lifestyle-related decline. In: Guccione AA, Wong RA, Avers D. *Geriatric Physical Therapy*. 3<sup>rd</sup> edition. Philadelphia, PA: Mosby; 27-37.
- Burkhead LM, Sapienza CM & Rosenbek JC. (2007). Laryngeal Elevation Achieved by Neuromuscular Stimulation at Rest. J Appl Physiol., 94:128-134.
- Campbell R, Polanksy R & McAdoo R. (1998). Deciphering Dysphagia with E-Stim. Advance for Directors in Rehabilitation, Nov: 7(11): 64-68.
- **Campbell R, Polanksy R & McAdoo R.** (1998). Promoting Laryngeal Elevation with E-Stim. Advance for Speech-Language Pathologists & Audiologists, June: 8(25) 24-27.
- Chae J, Bethoux F, Bohine T, Dobos L, Davis T & Friedl A. (1998). Neuromuscular Stimulation for Upper Extremity Motor and Functional Recovery in Acute Hemiplegia. *Stroke*, 29; 975-979.
- **Crary M.** (1995). A Direct Intervention Program for Chronic Neurogenic Dysphagia Secondary to Brainstem Stroke. *Dysphagia*, 10:6-18.
- Daniels SK, Brailey K, Priestly DH, Herrington LR, Weisberg LA & Foundas AL. (1998). Aspiration in patients with acute stroke. *Arch Phys Med Rehabil.*, 79:9-14.
- Dantas RO, Kern MK, Massey BT, Dodds WJ, Kahrilas PJ, Brasseur JG, Cook IJ & Lang IM. (1990). Effect of swallowed bolus variables on the oral and pharyngeal phases of swallowing. *Am J Physiol.*, 258:G675 G681.
- **Dorner, B.** (2002). It's tough to swallow: a practical approach to nutritional care of dysphagia. *Director,* Summer 10(3):107-10.
- **Doucet BM, Lam A & Griffin L.** (2012). Neuromuscular electrical stimulation for skeletal muscle function. Yale J Biol Med. Jun;85(2):201-215.
- Dumican MJ., Harper, K. (2022). ASHA Poster Presentation: Swallowing Outcomes in Clinically Based Patients Using Neuromuscular Electrical Stimulation (NMES)
- Elidan J, Gonen B, Shochina M & Gay I. (1990). Electromyography of the inferior constrictor and cricopharyngeal muscles during swallowing. *Ann Otol Rhinol Laryngol.*, 99:46-49.
- Enoka RM & Duchateau J. (2008) Muscle fatigue: what, why and how it influences muscle function. J Physiol. 586(Pt 1):11-23.
- Eriksson PO, Eriksson A, Ringqvist M & Thornell LE. (1982). Histochemical fibre composition of the human digastric muscle. *Arch Oral Biol.*, 27(3):207-15.
- Evans WJ. (2010) Skeletal muscle loss: cachexia, sarcopenia, and inactivity. Am J Clin Nutr, 91(4):1123S-1127S.

Evans WJ. (1999). Exercise training guidelines for the elderly. Med Sci Sports Exerc. 31(1):12-17.

- **Faghri PD, Rodgers MM, Glaser RM, Bors JG, Ho C & Akuthota P.** (1994). The Effects of Functional electrical stimulation on Shoulder Subluxation, Arm Function Recovery, and Shoulder Pain in Hemiplegic Stroke Patients. Arch Phys Rehabil. 75:73-79.
- **Fitzgerald GK, Piva SR & Irrgang JJ.** (2003). A Modified Neuromuscular Electrical Stimulation Protocol for Quadriceps Strength Training Following ACL Reconstruction. *JOSPT*, 33 (9): 492 501.
- Gargiulo P, Reynisson PJ, Helgason B, Kern H, Mayr W, Ingvarsson P, Helgason T & Carraro U. (2011). Muscle, tendons, and bone: structural changes during denervation and FES treatment. Neurol Res. Sep: 33(7):750-758.
- Groher ME. (1992). Dysphagia Diagnosis and Management. Boston: Butterworth-Heinemann.
- Henneman E, Clamann HP, Gillies JD & Skinner RD. (1974). Rank order of motorneurones within a pool, law of combination. J. Neurophysiology, 37:1338–1349.
- Henneman E, Somjen G & Carpenter DO. (1965). Functional significance of cell size in spinal motoneurons. J. Neurophysiology, 28:560–582.
- Hoffman MR, Mielens JD, Ciucci MR, Jones CA, Jiang JJ & McCulloch, TM. (2012). High resolution manometry of pharyngeal swallow pressure events associated with effortful swallow and the Mendelsohn maneuver. Dysphagia, 27(3), 418-426
- Hollinshead WH. (1982). Anatomy for Surgeons Volume 1: The Head and Neck. 3rd Ed. Philadelphia: Harper & Row, 367-371.
- Howden CW. (2004). Management of acid-related disorders in patients with dysphagia. American Journal of Medicine, Sept.6: 117 (5A): 44S-48S.
- **Hrycyshyn AW & Basmajian JV.** (1972). EMG of the oral stage of swallowing in man. *Am J Anatomy*, 133:333-340.
- Humbert I, Poletto CJ, Saxon KG, Kearney PR, Crujido L, Wright-Harp W, Payne J, Jeffries N, Sonies BC & Ludlow CL. (2006). The effect of surface electrical stimulation on hyo-laryngeal movement in normal individuals at rest and during swallowing. *J Appl Physiol.*, 101, 1657-1663.
- Jabre JF & Spellman NT. (1996). The demonstration of the size principle in humans using macro electromyography and precision decomposition. Muscle Nerve 19:338–341.
- Jacob P, Kahrilas P, Logemann JA, Shah V & Ha T. (1989). Upper esophageal sphincter opening and modulation during swallowing. *Gastroenterology*, 97:1469-1478.
- Kahrilas PJ, Dodds WJ, Dent J, Logemann JA & Shaker R. (1988). Upper esophageal sphincter function during deglutition. *Gastroenterology*, 95:52-62.
- Kibler WB. (1991). Role of the Scapula in the Overhead Throwing Motion. Contemp Orthop., 22:525.
- Kim J & Choi JY. (2016). The effect of subthreshold continuous electrical stimulation on the facial function of patients with Bell's palsy, Acta Oto-Laryngologica, 136:1, 100-105, DOI 10.3109/00016489.2015.1083121.
- Kleim J.A & Jones TA. (2008). Principles of experience-dependent neural plasticity: Implications for rehabilitation after brain damage. *Journal of Speech, Language, and Hearing Research: JSLHR*, 51(1), S225– S239.
- Korfage JA, Schueler YT, Brugman P & Van Eijden TM. (2001). Differences in myosin heavy-chain composition between human jaw-closing muscles and supra- and infrahyoid muscles. Arch Oral Biol. Sep: 46(9): 821-7.
- Lang IM & Shaker R. (1997). Anatomy and physiology of the upper esophageal sphincter. *Am J Med.*, 103:50S-55S.
- Langmore SE & McCulloch TM. (1997). Deglutition and its Disorders, San Diego: Singular Publishing.
- Langmore SE. (1991). Managing the Complications of Aspiration in Dysphagic Adults. New York: Thieme Medical Publishers.
- Llewellyn ME, Thompson KR, Deisseroth K, & Delp SL. (2010). Orderly recruitment of motor units under optical control in vivo. Nature Medicine, 16(10), 1161-1165.
- Lazarus, C, Logemann, JA, & Gibbons, P. (1993). Effects of maneuvers on swallowing function in a dysphagic oral cancer patient. *Head Neck*, 15(5), 419-424
- **Logemann JA.** (1999). Behavioral Management for Oropharyngeal Dysphagia. Folia Phoniatrica et Logopaedica, 51:199-212.
- Logemann JA. (1995). Dysphagia: Evaluation and Treatment. Folia Phoniatrica et Logopaedica, 47:140-164.
- Logemann JA. (1983). Evaluation and Treatment of Swallowing Disorders. San Diego: College-Hill Press.
- Ludlow CL, Humbert I, Saxon KG, Poletto CJ, Sonies BC & Crujido L. (2007). Effects of surface electrical stimulation both at rest and during swallowing in chronic pharyngeal dysphagia. *Dysphagia*, 22, 1-10.

- Ludlow CL, Wright-Harp W, Martin S & Safi M. (2010). Evaluation of submental stimulation in healthy volunteers. *Aerodigestive Tract Functions Conference*, Madison, WI.
- Lundy DS, Smith C, Colangelo L, Sullivan PA, Logemann JA, Lazarus CL, Newman LA, Murry T, Lombard L & Gaziano J. (1999). Aspiration: cause and implications. *Otolaryngol Head Neck Surg.*, 120: 474–478.
- Luschei ES1, Ramig LO, Baker KL & Smith ME. (1999) Discharge characteristics of laryngeal single motor units during phonation in young and older adults and in persons with Parkinson disease. J Neurophysiol. May;81(5): 2131-2139.
- Mann G, Hankey GJ & Cameron D. (2000). Swallowing disorders following acute stroke: prevalence and diagnostic accuracy. *Cerebrovasc Dis.*, 10: 380-386.
- Marchand, S. (2014). Coding and Billing for Therapy and Rehab.Sherry Marchand and Cross Country Education, LLC, 58-62.
- **McConnel FMS.** (1988). Analysis of pressure generation and bolus transit during pharyngeal swallowing. *Laryngoscope*, 98:71-78.
- McCullough, G, Kamarunas E, Mann, G, Schmidley J, Robbins J & Crary M. (2012). Effects of Mendelsohn Maneuver on Measures of Swallowing Duration Post-Stroke. *Top Stroke Rehabil.* 19(3): 234-243.
- **McCullough, G,** (2014). One step back and two steps up and forward: The superior movements of research defining the utility of the Mendelsohn Maneuver for improving UES function. *SIG 13 Perspectives on Swallowing and Swallowing Disorders*. 23:5-10.
- Medda BK, Lang IM, Layman R, Hogan WJ, Dodds WJ & Shaker R. (1994). Characterization and quantification of a pharyngo-UES contractile reflex in cats. *Am J Physiol.*, 167:(Gastrointest Liver Physiol 30)G972-G983.

**Medical Billing Terms.** Retrieved from https://www.all-things-medical-billing.com/medical-billing-codes.html **Miller AJ.** (1994). *Neuroscientific principles of swallowing.* Singular Publishing, London.

- Milner-Brown HS, Stein RB & Yemm R. (1973). The orderly recruitment of human motor units during voluntary isometric contractions. *J Physiol* 230:359–370.
- Molfenter SM & Steele CM (2014). Use of an anatomical scalar to control for sex-based size difference in measures of hyoid excursions during swallowing. Journal of Speech-Language Hearing Research, 57(3), 768-778. DOI: 10.1044/2014\_JSLHR-S-13-0152.

Morgan A, Ward E, Murdoch B, Kennedy B, et al. (2003). Incidence, characteristics, and predictive factors for dysphagia after pediatric traumatic brain injury. *Journal of HeadTrauma & Rehabilitation*, 18 (3): 239-251.

Mosler BJ, Jobe FW, Pink M, Perry J & Tibone J. (1992). EMG Analysis of the Scapular Muscles During a Shoulder Rehabilitation Program. *AM J Sports Med.*, 20(2):128.

Mu L & Sanders I. (1998). Neuromuscular organization of the human upper esophageal sphincter. Ann Otol Rhinol Larygnol., May: 107:370-7.

Mu L & Sanders I. (1996). The innervation of the human upper esophageal sphincter. Dysphagia, 11:234-238.

Mueller MJ & Maluf KS. (2002). Tissue Adaptation to Physical Stress: A Proposed "Physical Stress Theory" to guide physical therapist practice, education, and research. *Physical Therapy*, 82(4): 383 – 403.

- Murakami Y, Fukuda H & Kirchner JA. (1973). The cricopharyngeus muscle. An electrophysiological and neuropharmacological study. *Acta Otolaryngol.*, 311(Suppl):4-19.
- Nekl CG, Lintzenich CR, Leng X, Lever T, & Butler SG. (2012). Effects of effortful swallow on esophageal function in healthy adults. *Neurogastroenterology & Motility*, 24(3): 252-e108.
- Nguyen NP, Frank C, Moltz CC, Vos P, Smith HJ, Bhamidipati PV, Karlsson U, Nguyen PD, Alfieri A, Nguyen LM, Lemanski C, Chan W, Rose S, & Sallah S. (2006). Aspiration rate following chemoradiation for head and neck cancer: an underreported occurrence. *Radiother Oncol.*, 80: 302-306.
- Nguyen NP, Moltz CC, Frank C, Vos P, Smith HJ, Karlsson U, Nguyen LM, Rose S, Dutta S, Nguyen N & Sallah S. (2008). Long-term aspiration following treatment for head and neck cancer. *Oncology*, 74: 25-30.
- **Ogura et al. (**2022). Immediate Effects of Electrical Stimulation on Oropharyngeal Structure and Laryngeal Vestibular Closure: A Pilot Study in Healthy Subjects. +. Vol. 7, 20220033
- Olsson R, Nilsson H & Ekberg O. (1995). Simultaneous videoradiography and pharyngeal solid state manometry (videomanometry) in 25 nondysphagic volunteers. *Dysphagia*, 10:36-41.
- **Osborne SL.** (1951). The retardation of atrophy in man by electrical stimulation in muscles. *Archives of Physical Medicine*, 32:523–528.
- **Page EB.** (2006). Understanding the Modified Barium Swallow Study. *Radiographic Imaging CEU Source Manual.*
- Palmer J, Rudin N, Lara G & Crompton A. (1992). Coordination of mastication and swallowing. *Dysphagia*, 7:187-200.

- Park T, Kim Y, Ko DH & McCullough G. (2010) Initiation and Duration of Laryngeal Closure During the Pharyngeal Swallow in Post-Stroke Patients. *Dysphagia*, 25:177-182.
- **Perlman AL, Grayhack JP & Booth BM.** (1992). The relationship of vallecular residue to oral involvement reduced hyoid elevation and epiglottic function. *Journal of Speech and Hearing Research*, 35:734-741.
- Perlman AL, Palmer PM, McCulloch TM & VanDaele DJ. (1999). Electromyographic Activity from Human Laryngeal, Pharyngeal, and Submental Muscles During Swallowing. J. Appl. Physiol., 86:1663-1669.

Perlman AL. (1995). Quantification of certain parameters associated with the risk of aspiration. ASHA, 36.

Perlman AL & Schulze-Delrieu K. (1997). Deglutition and its Disorders, San Diego: Singular Publishing.

**Power ML, Hamdy S, Goulermas JY, Tyrrell PJ, Turnbull I & Thompson DG.** (2009) Predicting Aspiration After Hemispheric Stroke from Timing Measures of Oropharyngeal Bolus Flow and Laryngeal Closure. *Dysphagia*, 24:257-264.

Prentice W. (1994). Therapeutic Modalities in Sports Medicine, Oxford: Blackwell Scientific.

- **Ren J, Shaker R, Lang I & Sui Z.** (1995). Effect of volume, temperature and anesthesia on the pharyngo-UES contractile reflex in humans. *Gastroenterology*, 108:A677.
- Robbins JA, Gangnon R, Theis S, Kays S, Hewitt, A & Hind, J. (2005) The effects of lingual exercises on swallowing in older adults. *Journal of the American Geriatrics Society.*
- Rosenbek JC &Troche MS. (2013) Progressive Neurologic Disease and Dysphagia (Including Parkinson's Disease, Multiple Sclerosis, Amyotrophic Lateral Sclerosis, Myasthenia Gravis, Post-Polio Syndrome). In: Shaker R., Belafsky P., Postma G., Easterling C. (eds) Principles of Deglutition. Springer, New York, NY. https://doi.org/10.1007/978-1-4614-3794-9\_27.
- Selley WG, Ellis RE, Flack FC, Bayliss CR & Pearce VR. (1994). The synchronization of respiration and swallow sounds with videofluoroscopy during swallowing. *Dysphagia*, 9:162-167.
- Shaker R, Kern M, Bardan E, et al. (1997). Augmentation of deglutitive upper esophageal sphincter opening in the elderly by exercise. *Am J Physiol*, 272(35): G1518-O1522.
- Shaker R, Easterling C, Kern M, et al. (2002). Rehabilitation of swallowing by exercise in tube-fed patients with pharyngeal dysphagia secondary to abnormal UES opening. *Gastroenterology*, 122:1314-1321.
- **Shumway-Cook A & Woollacott MH.** (1995). Motor Control: Theory and Practical Applications. Baltimore: Williams and Wilkins.
- **Singh MA.** (2002). Exercise comes of age: rationale and recommendations for a geriatric exercise prescription. *J Gerontol A Biol Sci Med Sci.*,57(5):M262–M282.
- Smith J, Wolkove N, Colacone A & Kreisman H. (1989). Coordination of eating, drinking and breathing in adults. *Chest*, 96:578-582.
- Snyder-Macker L, Delitto A, Bailey S & Stralka SW. (1995). Quadriceps Femoris Muscle Strength and Functional Recovery After Anterior Cruciate Ligament Reconstruction: A Prospective Randomized Clinical Trial of Electrical Stimulation. J Bone Joint Surg., 77-A: 1166-1173.
- **Sproson L, Pownall S, Enderby P & Freeman J**. (2018). Combined electrical stimulation and exercise for swallow rehabilitation post-stroke: a pilot randomized control trial. *International Journal of Language & Communication Disorders*, 53(2): 405-417.
- Sundgren P, Maly P & Gullberg B. (1993). Elevation of the larynx on normal and abnormal cineradiogram. Br J Radiol., 66: 768–772.
- Thein LA & Greenfield BH. (1997). Impingement Syndrome and Impingement-Related Instability. In: Donatelli RA, editor. Physical Therapy of the Shoulder. 3rd Ed. New York: Churchill Livingstone Inc.
- Vivodtzev I, Debigaré R, Gagnon P, Mainguy V, Saey D, Dubé A, Paré MÈ, Bélanger M & Maltais F. (2012). Functional and muscular effects of neuromuscular electrical stimulation in patients with severe COPD: a randomized clinical trial. *Chest.* Mar; 141(3):716-725.
- **Watts CR.** (2013). Measurement of Hyolaryngeal Muscle Activation Using Surface Electromyography for Comparison of Two Rehabilitation Dysphagia Exercises. *Arch of Phys. Med. and Rehab*, 94:2542-8.
- Watts CR, Dumican MJ. (2018). The effect of transcutaneous neuromuscular electrical stimulation on laryngeal vestibule closure timing in swallowing. *BMC Ear, Nose and Throat Disorders*. 18:5.
- Wheeler-Hegland KM, Rosenbek JC, Sapienza CM. (2008). Submental sEMG and Hyoid Movement During Mendelsohn Maneuver, Effortful Swallow, and Expiratory Muscle Strength Training. *Journal of Speech, Language, Hearing Research*, 51:1072-1087.
- Williams PL & Warwick R editors. (1980). Grays Anatomy 36th British Ed. Philadelphia : WB Saunders Company.
- Yokoyama M, Mitomi N, Tetsuka N & Niimi S. (2000). Role of laryngeal movement and effect of aging on swallowing pressure in the pharynx and upper esophageal sphincter. *Laryngoscope*,110:434-439.

#### **Orders:**

- Should indicate the type of modality (ESP<sup>™</sup>)
- The area of the body where the modality will be used
- The reason/purpose for the treatment and the frequency and duration

Example: Speech and Language Pathology (SLP) to add to Dysphagia Plan of Care (POC) Effective Swallowing Protocol (ESP) to the suprahyoidal and facial musculature to increase oropharyngeal contraction and improve facial muscle tone 5 times a week for 30 days.

# Goal writing should address impairments that limit function and prevent complications. Goals should be:

- Measurable
  - Percentages are not always necessary, but when used indicate the number of stimuli, i.e., 80% accuracy (8/10)
- Lead to functional outcomes
  - Increased independence for ADL's
  - Less burden on caregiver and less cost to insurer
- Related to functional skills
- Patient specific customized to each patient's prior level of function

#### Sample Long-Term Goals:

Pt will be able to tolerate PO intake of choice without choking or aspirating.

Pt to display a safe and efficient swallow of (insert specific diet consistency) with/without compensatory strategies to meet nutrition and hydration needs.

#### Sample Short-Term goals:

(References) 1 tsp (teaspoon) = 5 ml or 5 cc 1ml = 1 cc 3 tsp (teaspoons)/1 tablespoon = 15 ml or 15 cc

Pt to perform \_\_\_\_\_ (e.g. 48) out of 60 effortful swallowing exercises using ESP to suprahyoidal musculature to improve hyolaryngeal excursion and secretion management, increase UES opening, improve laryngeal vestibule closure, reduce pharyngeal residue and decrease risk of penetration and aspiration during oral intake.

Daily note idea: Depending on the patient you might want to modify the exercises per 10-minute intervals based on the repetitions from the duty cycle (e.g. 5 seconds on/ 25 seconds off = 60 reps):

Ex: Pt performed 17 out of 20 Mendelsohn Maneuvers.

Pt performed 16 out of 20 Chin to chest with RPD with an effortful swallow.

Pt performed 15 out of 20 Chin to chest with RPD, jaw open/closed with an effortful swallow.

-Pt. achieved 48 out of 60 exercises, will progress to a 5/20 duty cycle next treatment due to patient achieving 80% accuracy with effortful swallows.

Increase accuracy of the Mendelsohn Maneuver from 50% (5/10) to 70% (7/10) to improve airway protection.

Pt will improve hyolaryngeal excursion by performing hyolaryngeal exercises [(Mendelsohn, Shaker, Masaka, Falsetto, etc. (insert percentage 80%, accuracy 4/5)] in combination with the SLP using ESP to the suprahyoidal musculature to facilitate muscle reeducation and strength to establish a more effective swallow.

Pt will tolerate \_\_\_\_\_ (15 or 30) minutes of ESP to suprahyoidal musculature to improve (hyolaryngeal excursion, manage secretions) and increase UES elongation and laryngeal excursion to reduce pyriform retention and decrease the risk of penetration and aspiration during oral intake.

Resident will be able to perform Mendelsohn Maneuver on command 80% (4/5) of the time in 5-10 minutes, with or without PO, to facilitate hyolaryngeal excursion, increase UES elongation and laryngeal vestibule closure to reduce pyriform retention and decrease the risk of penetration and aspiration during oral intake.

Pt will demonstrate the ability to manage oral secretions better by initiating cued volitional swallows every 2 minutes (15 swallows per session) during treatment of dysphagia with clear vocal quality.

Pt. will demonstrate the ability of AP propulsion and tongue base retraction with 1-2 tsp bolus (liquids, puree) during (meal, trial feedings, 4/5 attempts or bites).

Pt will increase initiation of swallow with presentation of 15cc/tablespoon size bolus from \_\_\_\_\_ (5) seconds to \_\_\_\_\_ (3) seconds or less to (reduce pharyngeal delay, reduce posterior bolus leakage) decreasing the risk of aspiration.

Pt will demonstrate a (therapeutic technique) for safe and efficient swallow of 15 cc/tablespoon or cup sip of \_\_\_\_\_ (thin) liquid with 90% effectiveness (9/10). (Therapeutic technique could include: Supraglottic swallow, Mendelsohn)

# May need to work with other disciplines (OT/PT) to assess/modify/adapt Restorative Postural Device in order to:

Achieve appropriate neck posture by (<u>degrees extension</u>) to promote (increased/adequate) pharyngeal space for functional epiglottic closure to decrease the risk of penetration and aspiration.

Facilitate upright head control by modifying a restorative posture device to decrease signs and symptoms of aspiration.

Pt will achieve appropriate head and neck posture to manage their secretions (anterior bolus loss) by wearing a restorative posture device for \_\_\_\_ (2) hours while sitting in chair or participating in direct or indirect dysphagia treatment.

Pt and/or caregiver to be instructed in restorative posture device application and be able to give correct return demonstration \_\_\_\_\_ (100%) of the time.

## IDDSI (International Dysphagia Diet Standardisation Initiative)<sup>1</sup>



## 8-Point Penetration-Aspiration Scale<sup>2</sup>

<u>Score</u>	Description of Events
1	Material does not enter airway.
2	Material enters the airway, remains above the vocal folds, and is ejected from the airway.
3	Material enters the airway, remains above the vocal folds, and is not ejected from the airway.
4	Material enters the airway, contacts the vocal folds, and is ejected from the airway.
5	Material enters the airway, contacts the vocal folds, and is not ejected from the airway.
6	Material enters the airway, passes below the vocal folds, and is ejected into the larynx or out of the airway.
7	Material enters the airway, passes below the vocal folds, and is not ejected from the trachea despite effort.
8	Material enters the airway, passes below the vocal folds, and no effort is made to eject.

<sup>&</sup>lt;sup>1</sup> © The International Dysphagia Diet Standardisation Initiative 2019 @ https://iddsi.org/framework. Licensed under the CreativeCommons Attribution Sharealike 4.0 License. Derivative works extending beyond language translation are NOT PERMITTED.

<sup>&</sup>lt;sup>2</sup> Rosenbek, JC, Robbins, J, Roecker EV, Coyle, JL, & Woods, JL. A Penetration-Aspiration Scale. <u>Dysphagia</u> 11:93-98, 1996.

## Functional Oral Intake Scale<sup>3</sup>

#### TUBE DEPENDENT (levels 1-3)

- 1 No oral intake
- 2 Tube dependent with minimal/inconsistent oral intake
- 3 Tube supplements with consistent oral intake

#### TOTAL ORAL INTAKE (levels 4-7)

- 4 Total oral intake of a single consistency
- 5 Total oral intake of multiple consistencies requiring special preparation
- 6 Total oral intake with no special preparation, but must avoid specific foods or liquid items
- 7 Total oral intake with no restrictions

## Waxman Dysphagia Severity Scale<sup>4</sup>

Rating	Explanation
0	Normal swallowing mechanism.
1	Minimal dysphagia—video swallow shows slight deviance from a normal swallow. Patient may report a change in sensation during swallow. No change in diet is required.
2	Mild dysphagia—oropharyngeal dysphagia present, which can be managed by specific swallow suggestions. Slight modification in consistency of diet may be indicated.
3	Mild-moderate dysphagia—potential for aspiration exists but is diminished by specific swallow techniques and a modified diet. Time for eating is significantly increased; thus supplemental nutrition may be indicated.
4	Moderate dysphagia—significant potential for aspiration exists. Trace aspiration of one or more consistencies may be seen under videofluoroscopy. Patient may eat certain consistencies by using specific techniques to minimize potential for aspiration and/or to facilitate swallowing. Supervision at mealtimes required. May require supplemental nutrition orally or via feeding tube.
5	Moderately severe dysphagia—patient aspirates 5% to 10% on one or more consistencies, with potential for aspiration on all consistencies. Potential for aspiration minimized by specific swallow instructions. Cough reflex absent or nonprotective. Alternative mode of feeding required to maintain patient's nutritional needs. If pulmonary status is compromised, "nothing by mouth" may be indicated.
6	Severe dysphagia—more than 10% aspiration for all consistencies. "Nothing by mouth"

<sup>&</sup>lt;sup>3</sup> Crary MA, Carnaby-Mann GD, Groher ME. Initial psychometric assessment of a functional oral intake scale for dysphagia in stroke patients. *Arch Phys Med Rehabil* 2005;86:1516-1520.

recommended.

<sup>&</sup>lt;sup>4</sup> Waxman, Michael J., et al. "Nutritional aspects and swallowing function of patients with Parkinson's disease." *Nutrition in clinical practice* 5.5 (1990): 196-199.

System	Synkinesis	Rate the degree of INVOLUNTARY MUSCLE CONTRACTION associated with each expression	Line of the second sec	0 1 2 3	0 1 2 3	0 1 2 3	0 1 2 3	0 1 2 3		Synkinesis score: Total	Composite score
nybrook Facial Grading	Symmetry of Voluntary Movement	Degree of muscle EXCURSION compared to normal side	Standard Kovement entitiete Standard Expressions Expr	Forehead 1 2 3 4 5	Gentle eye closure (OCS) 1 2 3 4 5	Open mouth 1 2 3 4 5	Snarl 1 2 3 4 5	Lip Pucker 1 2 3 4 5	Local Lo	Voluntary movement score: Total × 4	Vol Resting O Synk Score Synk Score
Sun	<b>Resting Symmetry</b>	Compared to normal side	Eye (choose one only) normal 0 narrow 1 wide 1 eyelid surgery 1 eyelid surgery 1 cheek (naso-lablal fild) normal 0	less pronounced 1 more manunced 1		Mouth normal D	corner dropped corner pulled up/out 1		Resting symmetry Total × 5	Patient's name	Dx Date 1

House-Brackmann Facial Nerve Grading System <sup>2</sup>
<b>Grade I - Normal</b> Normal facial function in all areas
<b>Grade II - Slight Dysfunction</b> Gross: slight weakness noticeable on close inspection; may have very slight synkinesis At rest: normal symmetry and tone Motion: <u>forehead</u> - moderate to good function; <u>eye</u> - complete closure with minimum effort; <u>mouth</u> - slight asymmetry.
Grade III - Moderate Dysfunction Gross: obvious but not disfiguring difference between two sides; noticeable but not severe synkinesis, contracture, and/or hemi-facial spasm. At rest: normal symmetry and tone Motion: <u>forehead</u> - slight to moderate movement; <u>eve</u> - complete closure with effort; <u>mouth</u> - slightly weak with maximum effort.
<b>Grade IV - Moderate Severe Dysfunction</b> Gross: obvious weakness and/or disfiguring asymmetry At rest: normal symmetry and tone Motion: <u>forehead</u> - none; <u>eve</u> - incomplete closure; <u>mouth</u> - asymmetric with maximum effort.
<b>Grade V - Severe Dysfunction</b> Gross: only barely perceptible motion At rest: asymmetry Motion: <u>forehead</u> - none; <u>eve</u> - incomplete closure; <u>mouth</u> - slight movement
<b>Grade VI - Total Paralysis</b> No movement
<sup>1</sup> Sunnybrook Health Sciences Centre (2019). https://sunnybrook.ca/uploads/FacialGradingSystem.pdf <sup>2</sup> House JW, Brackmann DE (1985). Facial nerve grading system. Otolaryngol. Head Neck Surg, [93] 146–147.

## EATING ASSESSMENT TOOL (EAT-10)

Name: \_\_\_\_\_ EMR# \_\_\_\_\_

Birthdate/Age: \_\_\_\_\_ Todays Date: \_\_\_\_\_

The purpose of the EAT-10 questions is to help measure swallowing problems. Answer each question by circling the number that matches how bad you feel the problem is for you.

To what degree to you experience the							
following problems?	0 = No problem 4 = Severe problem						
Circle an answer between 0 and 4							
1. My swallowing problem has caused me to	0	1	2	3	4		
lose weight.							
2. My swallowing problem interferes with my	0	1	2	3	4		
ability to go out for meals.							
3. Swallowing liquids takes extra effort.	0	1	2	3	4		
4. Swallowing solids takes extra effort.	0	1	2	3	4		
5. Swallowing pills takes extra effort.	0	1	2	3	4		
6. Swallowing is painful.	0	1	2	3	4		
7. The pleasure of eating is affected by my	0	1	2	3	4		
swallowing.							
8. When I swallow food sticks in my throat.	0	1	2	3	4		
9. I cough when I eat.	0	1	2	3	4		
10.Swallowing is stressful	0	1	2	3	4		
		TOTAL	FAT 40	<b>C</b>			

Add up the sum of the numbers you circled for a TOTAL EAT-10 Score:

If your score is greater than 3 you may have swallowing problems. We suggest that you share your EAT-10 results with your doctor.

Reference: Belafky PC, Mouadeb DA, Rees CJ, Pryor JC, Postma GN, Alen J and Leonard RJ. Validity and reliability of the Eating Assessment Tool (EAT-10). Ann Otol Rhinol Laryngol 2008; 117(12):919-924.



#### Today's Date: \_\_\_\_\_ Certification Number: \_\_\_\_\_

## ESP™ Dysphagia Treatment System Order Form

Email your order to:info@ampcarellc.comFax: 817-348-8830Phone: 682-561-2444Orders online: www.ampcarellc.com/orders

#### Bill to:

#### Ship to (if different than Bill to information):

Customer/Facility Name	Customer/Facility Name
Accounts Payable Contact	Attention
Address	A deleses
Address	Address
City State 7in	City State Zin
Phone Fax	Phone
Email	Email

#### **Method of Payment**

	aymene		
Invoice		PO#	Authorized Signature
Visa	MasterCard	Account#	Exp. Date
AmEx	Discover	Cardholder Name	Security Code
		Signature	Billing Zip Code

#### **Product Selection**

	Product	roduct			
Qty	Number	Product Description	Price	Total	
AMPC	ARE ESP THERA	PY SYSTEM			
		ESP™ Therapy System Kit {Ampcare ES™ Unit, user manual, batteries, lead wires and tester,			
	50706LT	10 pkg Large/Adult E Series Electrodes, Restorative Posture Device (RPD) and replacement pad	\$ 799.00		
	50706ST	\$ 799.00			
		set, self-adhering wrap, shoulder-strap carrying bag}			
LARG	E (Adult) ELECTF	RODES 1.5" x 1.75":			
	50709LT-1	10 Pack E Series Large Electrodes (4/pkg, 10 pkg/pouch) – 1 pouch	\$ 139.50		
	50709LT-3	30 Pack E Series Large Electrodes (4/pkg, 10 pkg/pouch) – 3 pouches	\$ 388.50		
	50709LT-5	50 Pack E Series Large Electrodes (4/pkg, 10 pkg/pouch) – 5 pouches	\$ 597.50		
	\$ 995.00				
SMAL	L (Youth) ELECT	RODES 1.0" x 1.25":	.=		
	50709ST-1	10 Pack E Series Small Electrodes (4/pkg, 10 pkg/pouch) – 1 pouch	\$ 139.50		
	50709ST-3	30 Pack E Series Small Electrodes (4/pkg, 10 pkg/pouch) – 3 pouches	\$ 388.50		
	50709ST-5	50 Pack E Series Small Electrodes (4/pkg, 10 pkg/pouch) – 5 pouches	\$ 597.50		
	50709ST-10	100 Pack E Series Small Electrodes (4/pkg, 10 pkg/pouch) – 10 pouches	\$ 995.00		
ACCES	SORIES				
[	50707	RPD (Restorative Posture Device)	\$ 100.00		
	50708	RPD Replacement Pad Set	\$ 25.00		
	F0710LT	Ampcare ES™ (ES unit, user manual, batteries, lead wires, lead wire tester,	ć 740.00		
	50710L1	10 pkg Large/Adult E Series Electrodes, and Ampcare shoulder-strap carrying bag)	\$ 749.00		
	Ampcare ES <sup>™</sup> (ES unit, user manual, batteries, lead wires, lead wire tester,				
	50710ST 10 pkg <u>Small/</u> Youth E Series Electrodes, and Ampcare shoulder-strap carrying bag)				
[	50712	Replacement Lead Wire Set (2) for Ampcare ES™ with Lead Wire Tester	\$ 40.00		
	50713	Lead Wire Tester	\$ 15.00		
	50714	Ampcare Shoulder Strap Carrying Bag	\$ 35.00		
Sales ta Taxes,	ax applies where r shipping and hand	equired by law. If you are tax exempt, please fax your tax exemption certificate along with your order. dling will be added to the final invoice.	Subtotal		

## Ampcare's ESP<sup>™</sup> Placements

## 1. Submental



- Lead wires will extend beyond the front of the mandible
- Electrodes should <u>not</u> touch one another or the mandible
- Do not place lower than the thyroid notch

#### **General Rules:**

## 2. Facial



- Back electrode should be no higher than the acoustic meatus to prevent eye twitch
- Front electrode can cover both buccinator and orbicularis oris
- Consider smaller electrodes if you want to target specific muscles

Clean the area(s) with soap and water or alcohol wipe and let dry (approx. 45 seconds) Electrodes should never touch one another Place the black lead over the weaker side or musculature Use Posture Device or coban strap to keep electrodes in place

#### Sample Verbal Instructions for the ESP:

I would like to use this device to rehabilitate some of your swallowing muscles, would that be alright?

I will place (two or four) pads on your skin (either under your chin and or face) and connect them to the device with this lead wire. The pads may feel a little cold when it first touches your skin.

I will place this posture device under your chin, on your chest and fasten around the back of your neck or use this coban strap around your chin.

As I am slowly turning up the vibrations please <u>tell me when you feel a tingle</u> under the pads.

<u>Take as much</u> "tingling" <u>as you can</u> in attempt to make the muscle tighten and wake up. "<u>The</u> <u>more</u> tingling <u>the better</u>."

I will ask you to perform exercises while the vibration is on.

#### Preset Programming ESP 1 (50 phase duration) or ESP 2 (250 phase duration):

If you are using a Preset program you will only be able to modify the intensity and the cycle "off" time.

## The Keypad



#### To program a treatment session:

- 1. Press the On/Off button to turn on the power to the unit.
- 2. Hold the program button Prog down for 3 seconds until either ESP 1 or ESP 2 is shown on the screen.
- 3. If you do not want to modify this program you may proceed to step 7.
- 4. To modify the cycle off time press the "Menu up" or "Menu down" button once. The off time will begin blinking to allow it to be modified.
- 5. Press the minus "-" or plus "+" button until the recommended cycle off time is selected. This can be changed to 25, 20 or 15 seconds to lessen the amount of rest time between stimulations and therefore increase the amount of exercise repetitions.
- 6. Once the parameter has been selected press the enter button

Enter to accept it.

#### To start a treatment session:

- 7. Use the up and down arrows until the prescribed intensity is visible in the LCD (left arrows for channel 1 and right arrows for channel 2).
- 8. The timer will begin counting down automatically once the intensity is increased and the LED Light Indicator will blink green when the unit is administering current.
- 9. Once the treatment time has reached zero the unit will automatically power off.
- 10. Once ESP 1 or 2 has been programmed, this will be the default until it is reprogrammed.

#### Remember the mechanical causes/impairments that you can address with ESP:

- 1. Impaired lip closure
- 2. Poor mastication and bolus control
- 3. Increased oral pocketing/retention
- 4. Reduced tongue base retraction
- 5. Delayed timing and initiation of a swallow
- 6. Decreased hyolaryngeal excursion

- 7. Limited epiglottis retroversion
- 8. Delayed/incomplete laryngeal vestibule closure
- 9. Retention in the vallecular/pyriform sinus
- 10. Decreased pharyngeal contraction
- 11. Reduced opening of the PES/UES
- 12. Impaired conduction of neural pathways











### Cranial Nerves Involved in Swallowing and Their Evaluation Techniques

The term swallow reflex is deceptive in the sense that it is an elicited response requiring both sensory and motor information to initiate the swallow response.

<u>CN V (5) – Trigeminal</u> – conveys sensation to the face and motor to the soft palate (tensor veli palatini), pharynx and muscles of mastication. Primary interests to SLPs are Maxillary and Mandibular branches. (Oral and Pharyngeal Phase)

How to test:

- 1. Test with a wisp of cotton gently to the nostrils which should result in the wrinkling of the nose.
- 2. Bulk of masseter and temporalis muscle is tested by palpation of the muscles while the patient clenches their teeth. (Oral)
- 3. After bulk is determined, strength of muscles should be tested by attempting to open the jaw. Jaw will deviate to the weak side. (Oral)
- 4. Look inside mouth, weakness of tensor veli palatini may manifest itself as tilting of the uvula to the weak side. (Oral)
- 5. Palpate the hyoid while asking the patient to swallow and assess the amount of hyolaryngeal elevation. {Pharyngeal}

<u>CN VII (7) – Facial</u> – conveys sensation (taste) to the anterior 2/3 of the tongue, soft palate and motor to the facial muscles. (Oral Phase)

Taste is tested using substances that are sweet (sugar), sour (lemon swab) or salty (salt).

Test by inspection of facial symmetry and then individual muscle for strength.

"Please raise your eyebrows" (frontalis).

"Please close your eyes tightly and stop me from opening them" (orbicularis oculi).

"Pretend to blow out some candles" (orbicularis oris).

"Puff out your cheeks" (buccinator). Try to push the air out – while keeping your lips sealed. (orbicularis oris).

<u>CN IX (9) – Glossopharyngeal</u> – conveys sensation to the posterior 1/3 of the tongue, soft palate and pharynx and motor to the pharyngeal musculature. (Oral and Pharyngeal Phase)

Testing this musculature is accomplished by testing the gag reflex. The gag reflex is a true reflex. Usually this is triggered by placing a tongue depressor at the back of the tongue or along the soft palate. Your gag response comes from the brainstem function which is why this reflex is often affected with brainstem strokes. Using a lemon swab to the back of the tongue would allow you to test for both sensory and motor. (Oral and Pharyngeal)

<u>CN X (10) – Vagus</u> – conveys sensation to the trachea, larynx and pharynx and motor to the soft palate (levator veli palatini), larynx and pharynx along with other areas. {Pharyngeal Phase}

- 1. Test by inspecting palate, it will be lower and less arched. Say "ahhh" and you will see deviation to the normal side.
- 2. Test posterior pharyngeal wall gag reflex
- 3. Bilateral weakness may sound breathy or hyper nasal.
- 4. Unilateral may have a hoarse voice.

This nerve also influences esophageal motility and UES closure and opening.



<u>CN XI (11) – Accessory</u> – conveys motor to the shoulder, neck and soft palate (musculus uvulae). {Pharyngeal Phase}

Inspect for asymmetry and assess power.

"Please shrug your shoulders up and stop me from pushing them down."

"Please turn your head to the right and stop me from pushing it back." Feel the left sternocleidomastoid at the same time. Repeat on the other side.

<u>CN XII (12) – Hypoglossal</u> – conveys motor to the tongue muscles. (Oral Phase)

Test tongue motor function by protruding tongue, push tongue against cheek. The tongue will deviate to the side of the lesion.

Check for tongue fasciculations - they can be the first clue to impaired neurological function.

One muscle(s) and/or one nerve involvement can affect the normal swallow function. Assessment: Pt. appears to present with possible impairment of \_\_\_CN due in part to deficit in \_\_\_ test.

	Orar i flase	
Swallowing Problem	Muscle Group	Cranial Nerve
Acceptance	Orbicularis Oris, Zygomaticus	CN VII
Anterior loss	Orbicularis Oris	CN VII
Mastication	Masseter, Pterygoids, and Temporalis	CN V
Oral pocketing or scattering	Buccinator and intrinsic and extrinsic lingual muscles	CN VII, XII
Anterior-posterior propulsion	Intrinsic and extrinsic lingual muscles – Genioglossus, Styloglossus	CN XII
Poor lingual pressure to drive the bolus through the PES/UES	Intrinsic and extrinsic lingual muscles – Hyoglossus	CN XII

**Oral Phase** 

Swallowing Problem	Muscle Group	Cranial Nerve
Poor velopharyngeal seal	Tensor Veli Palatini, Pharyngeal Constrictors, Levator Veli Palatini, Musculus Uvulae	CN V, IX, X, XI
Decreased closure of the larynx	Aryepiglottic, Lat. Cricoarytenoids, Transverse Arytenoid	CN X, XI
Weak pharyngeal constrictor contraction	Superior, Middle, and Inferior Pharyngeal Constrictors	CN IX, X, XI
Decreased anterior superior elevation of the hyolaryngeal complex	Anterior Digastric and Mylohyoid	CN V
Failure of opening the UES/PES	Inferior Pharyngeal Constrictor, Cricopharyngeus, Superior Longitudinal Esophageal muscle	CN IX, X, XI

#### International Classification of Diseases (ICD-10) – Dysphagia Diagnosis Codes

- R13.11 Dysphagia, oral phase
- R13.12 Dysphagia, oropharyngeal phase
- R13.13 Dysphagia, pharyngeal phase
- R13.14 Dysphagia, pharyngoesophageal phase